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Phenylalanin ELISA Kit

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1. INTENDED USE	16
2. PRINCIPLE OF THE TEST	16
3. MATERIAL SUPPLIED	17
4. MATERIAL REQUIRED BUT NOT SUPPLIED	17
5. PREPARATION AND STORAGE OF REAGENTS	18
6. PRECAUTIONS	19
7. SPECIMEN COLLECTION AND PREPARATION	19
8. ASSAY PROCEDURE	20
<i>SAMPLE PREPARATION PROCEDURE</i>	20
<i>TEST PROCEDURE</i>	21
9. EVALUATION OF RESULTS	22
<i>EXPECTED RESULTS</i>	23
10. PERFORMANCE CHARACTERISTICS	24
<i>CROSS REACTIVITY</i>	24
<i>PRECISION AND REPRODUCIBILITY</i>	24
<i>SENSITIVITY</i>	25
<i>RECOVERY</i>	25
<i>LINEARITY</i>	25
<i>CORRELATION WITH HPLC-MS</i>	26
11. LIMITATIONS	26
12. REFERENCES	26
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	27

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12. LITERATUR

- Badawy AA, Dougherty DM, Richard DM. Specificity of the Acute Tryptophan and Tyrosine Plus Phenylalanine Depletion and Loading Tests I. Review of Biochemical Aspects and Poor Specificity of Current Amino Acid Formulations. *Int J Tryptophan Res.* 2010 Jan 1;2010(3):23-34.
- Badawy AA, Dougherty DM, Richard DM. Specificity of the Acute Tryptophan and Tyrosine Plus Phenylalanine Depletion and Loading Tests Part II: Normalisation of the Tryptophan and the Tyrosine Plus Phenylalanine to Competing Amino Acid Ratios in a New Control Formulation. *Int J Tryptophan Res.* 2010;3:35-47.
- Robinson OJ, Standing HR, DeVito EE, Cools R, Sahakian BJ. Dopamine precursor depletion improves punishment prediction during reversal learning in healthy females but not males. *Psychopharmacology (Berl).* 2010 Aug;211(2):187-95. Epub 2010 May 22.

- Singhatiraj E, Ngamreungphong S, Nugent K. Acute dopamine depletion syndrome. J Am Geriatr Soc. 2010 Aug;58(8):1600-1. doi: 10.1111/j.1532-5415.2010.02987.x.
- van Spronsen FJ, de Groot MJ, Hoeksma M, Reijngoud DJ, van Rijn M. Large neutral amino acids in the treatment of PKU: from theory to practice. J Inherit Metab Dis. 2010 Oct 26. [Epub ahead of print].

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Manual

Phenylalanine ELISA Kit

*For the determination of phenylalanine in human EDTA plasma
and serum*

For research use only



Valid from 24.11.2010

K 7014

96

1. INTENDED USE

The Phenylalanine ELISA Kit is intended for the quantitative determination of phenylalanine in human EDTA plasma and serum. It is for research use only.

2. PRINCIPLE OF THE TEST

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for phenylalanine derivatization. Afterwards, the treated samples and the polyclonal phenylalanine antiserum are incubated in wells of a microtiter plate coated with a phenylalanine-derivative (tracer). During the incubation period, the target phenylalanine in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The phenylalanine in the sample displaces the antibodies out of the binding to the tracer. Therefore, the concentration of the tracer-bound antibody is inverse proportional to the phenylalanine concentration in the sample. During the second incubation step, a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-phenylalanine antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the phenylalanine concentration in the sample; this means, high phenylalanine concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Phenylalanine present in the patient samples is determined directly from this curve.

3. MATERIAL SUPPLIED

Catalog No	Content	Kit Components	Quantity
K7014MTP	PLATE	One holder with precoated strips	12 x 8 wells
K7014ST	STD	Standards (0, 25, 50, 100, 200, 400 μ M), diluted in reaction buffer (ready to use)	6 x 1 vial
K7014KO1 K7014KO2	CTRL 1 CTRL 2	Controls, diluted in reaction buffer (ready to use)	2 x 1 vial
K7014WP	WASHBUF	Wash buffer concentrate (10-fold)	2 x 100 ml
K7014AK	AB	Anti-phenylalanine antibody (lyophilized)	2 x 1 vial
K7014VR	ABBUF	Antibody dilution buffer concentrate (10-fold)	2 x 1 ml
K7014K	2.AB	POD antibody (concentrate)	120 μ l
K7014CSP	2.ABDIL	Conjugate stabilizing buffer	24 ml
K7014RP	DERBUF	Reaction buffer	2 x 25 ml
K7014DR	DER	Derivatization reagent	2 x 25 mg
K7014LM	DMSO	Dimethylsulfoxide (DMSO)	3.5 ml
K7014SL	CODIL	Dilution buffer for coupling	2 x 28 ml
K7014TMB	SUB	TMB substrate	25 ml
K7014AC	STOP	Stop solution	15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Double distilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 10000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm
(reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with aqua bidest. **1:10** before use (**100 ml WASHBUF + 900 ml aqua bidest.**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The WASHBUF is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- **Standards (STD) and controls (CTRL1, CTRL2) are already diluted** in reaction buffer (DERBUF). Store standards and controls frozen at -20°C, thaw before use in the test, and re-freeze immediately after use. Standards and controls can be re-frozen up to 3 times.
- **DMSO** could crystallize at 4°C. Dissolve the crystals at 20-25°C in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (25.0 mg) in 1.5 ml DMSO.** Put the vial on a horizontal shaker for 5 min. Discard any rest of the reagent after use. DER must be **prepared immediately before use.** The ELISA kit can be separated into two performances by providing two DER vials. Please note: **DMSO attacks all plastics but not polypropylene products and laboratory glass.**
- Dilute the **antibody dilution buffer concentrate (ABBUF)** with aqua bidest. **1:10** before use (**1 ml concentrate + 9 ml aqua bidest.**), mix well. The ELISA kit can be separated into two performances by providing 2 x 1 ml antibody dilution buffer concentrate. The diluted buffer is stable over a longer period. It can be stored at **2-8°C for 4 weeks.**
- Dissolve the **anti-phenylalanine antibody (AB) in 6 ml of antibody dilution buffer.** The ELISA kit can be separated into two performances by providing two AB vials. Diluted anti-phenylalanine antibody is stable over a longer period. It can be stored at **2-8°C for 4 weeks.**

- Dilute the **POD antibody (2.AB) 1:200** with conjugate stabilizing buffer (2.ABDIL), **e.g. 110 µl 2.AB + 22 ml 2.ABDIL, prepare only the required amount.** The undiluted POD antibody (2.AB) is stable at **2-8°C** until the expiry date stated on the label. Diluted POD antibody (2.AB) is not stable over a longer period. It can be stored at **2-8°C for only 5 days.**
- All other test reagents are ready for use. Test reagents are stable until the expiry date (see label of test package) when stored at 2-8°C.

6. PRECAUTIONS

- For research use only.
- Human materials used in kit components were tested and found to be negative for HIV and Hepatitis B. However, for safety reasons all kit components should be treated as if potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

7. SPECIMEN COLLECTION AND PREPARATION

EDTA plasma and serum

- Venous fasting blood is suited for this test system. Blood samples are stable for one week at 2-8°C. For longer storage blood samples should be frozen at -20°C.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- The EDTA plasma and serum samples are diluted for derivatization.
Samples with visible amounts of **precipitates** should be **centrifuged** at least for 5 min at 10000 x g. The resulting supernatant is used in the assay.
- For sample preparation, a derivatization reagent (DER) for coupling of phenylalanine is added (details are given in the sample preparation procedure).

8. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according to the enclosed manual.

Sample preparation procedure

Dilute **EDTA plasma and serum samples** with reaction buffer by **factor 1:100**, i.e. **10 µl sample + 990 µl reaction buffer (DERBUF)**.

Derivatization of standards (STD), controls (CTRL) and diluted samples (SAMPLE) is carried out in single analysis.

1.	Bring all reagents and samples to room temperature (18-26°C).
2.	Add 100 µl of ready to use standards (STD) , 100 µl of ready to use controls (CTRL) and 100 µl of diluted samples (SAMPLE) in the corresponding vial.
3.	Add 50 µl of freshly prepared derivatization reagent (DER) into each vial (standards, controls and samples), mix well and incubate for 45 min on a shaker (180-240 rpm) at room temperature (18-26°C) .
4.	Afterwards add 1000 µl of dilution buffer (CODIL) into each vial, mix well and incubate for 45 min on a shaker (180-240 rpm) at room temperature (18-26°C) .

2 x 100 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5.	Mark the positions of standards (STD)/ controls (CTRL)/ samples (SAMPLE) in duplicate on a protocol sheet.
6.	Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
7.	Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution.
8.	For the analysis in duplicate, take 2 x 100 µl of standards (STD) / controls (CTRL) / samples (SAMPLE) out of the vial and add into the respective wells of the microtiter plate (PLATE).
9.	Add 100 µl diluted anti-phenylalanine antibody (AB) into each well. Cover the plate tightly.
10.	Incubate overnight (15-20 hours) at 2-8°C .
11.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution.
12.	Add 200 µl diluted POD antibody (2. AB) into each well.
13.	Cover plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal shaker (180-240 rpm).
14.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution.
15.	Add 200 µl of TMB substrate (SUB) into each well.

- | |
|--|
| 16. Incubate for 8-14 min at room temperature (18-26°C) in the dark*. |
| 17. Add 100 µl of stop solution (STOP) into each well, mix thoroughly. |
| 18. Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used. |

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

9. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, samples, and with correct sample treatment), standards, controls, and samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results**.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).


Plausibility of the measured pairs of values should be examined before automatically evaluating the results. If this option is not available within the used program, the pairs of values should be controlled manually.

Controls

Control samples or plasma pools should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from the calibration curve. In the following, an example of a calibration curve is given. Do not use it for the calculation of your results.

Example of calibration curve

Expected results

Based on internal studies with plasma samples of evidently healthy persons (**n=40**) a mean value of 97 $\mu\text{mol/l}$ was calculated. The standard deviation was 19.7 $\mu\text{mol/l}$.

Plasma mean value $\pm$ 2 x standard deviation: **97 $\pm$ 39.4 $\mu\text{mol/l}$**

Normal range: **57.6 - 136.4 $\mu\text{mol/l}$**

We recommend each laboratory to develop its own normal range. The values mentioned above are only for orientation and can deviate from other published data.

10. PERFORMANCE CHARACTERISTICS

Cross reactivity

L-Thyrosine < 0.02 %

Precision and reproducibility

EDTA plasma:

Intra-assay (n=10)		
sample	phenylalanine [$\mu\text{mol/l}$]	coefficient of variation (CV) [%]
1	47.8	8.3
2	83.61	10.2

Inter-assay (n=6)		
sample	phenylalanine [$\mu\text{mol/l}$]	coefficient of variation (CV) [%]
1	45.5	10.9
2	79.47	5.1

Sensitivity

The sensitivity was set as $B_0 + 2SD$. The zero-standard was measured 48 times.

sample	phenylalanine mean value [OD]	2 x standard deviation (SD)	detection limit [$\mu\text{mol/l}$]
zero-standard	1.9	0.1	16.9

Recovery

One sample was spiked with different phenylalanine concentrations and measured using this assay. The analytical recovery rate was determined by the expected and measured phenylalanine levels. The expected levels were calculated as the sum of the measured phenylalanine concentration in the original sample and the spiked phenylalanine amount. The mean recovery rate for all concentrations was 103.07 % (n=10).

EDTA plasma:

spike [$\mu\text{mol/l}$]	phenylalanine measured [$\mu\text{mol/l}$]	phenylalanine expected [$\mu\text{mol/l}$]	recovery [%]
0	x = 40.9	x	100.0
20	62.0	40.9+20 = 60.9	101.8
50	97.7	40.9+50 = 90.9	107.5

Linearity


The linearity of the ELISA was determined by the dilution of a spiked patient sample. The mean linearity was 100.8 % (n=10).

EDTA-plasma

dilution	measured [$\mu\text{mol/l}$]	expected [$\mu\text{mol/l}$]	recovery [%]
original	97.7	x	100.0
2+1	62.0	65.1	95.1
1+1	52.4	48.8	107.4

Correlation with HPLC-MS

14 samples were measured with this ELISA and HPLC-MS. The correlation was $r = 0.918$.

11. LIMITATIONS

Hemolytic and lipemic samples may give erroneous results. Do not measure hemolytic and lipemic samples.

12. REFERENCES

- Badawy AA, Dougherty DM, Richard DM. Specificity of the Acute Tryptophan and Tyrosine Plus Phenylalanine Depletion and Loading Tests I. Review of Biochemical Aspects and Poor Specificity of Current Amino Acid Formulations. *Int J Tryptophan Res.* 2010 Jan 1;2010(3):23-34.
- Badawy AA, Dougherty DM, Richard DM. Specificity of the Acute Tryptophan and Tyrosine Plus Phenylalanine Depletion and Loading Tests Part II: Normalisation of the Tryptophan and the Tyrosine Plus Phenylalanine to Competing Amino Acid Ratios in a New Control Formulation. *Int J Tryptophan Res.* 2010;3:35-47.

- Robinson OJ, Standing HR, DeVito EE, Cools R, Sahakian BJ. Dopamine precursor depletion improves punishment prediction during reversal learning in healthy females but not males. *Psychopharmacology (Berl)*. 2010 Aug;211(2):187-95. Epub 2010 May 22.
- Singhatiraj E, Ngamreungphong S, Nugent K. Acute dopamine depletion syndrome. *J Am Geriatr Soc*. 2010 Aug;58(8):1600-1. doi: 10.1111/j.1532-5415.2010.02987.x.
- van Spronsen FJ, de Groot MJ, Hoeksma M, Reijngoud DJ, van Rijn M. Large neutral amino acids in the treatment of PKU: from theory to practice. *J Inherit Metab Dis*. 2010 Oct 26. [Epub ahead of print].

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for research use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number