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Tryptophan ELISA Kit

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1. INTENDED USE	18
2. INTRODUCTION	18
3. PRINCIPLE OF THE TEST	19
4. MATERIAL SUPPLIED	20
5. MATERIAL REQUIRED BUT NOT SUPPLIED	21
6. PREPARATION AND STORAGE OF REAGENTS	21
7. PRECAUTIONS	22
8. SPECIMEN COLLECTION AND PREPARATION	22
9. ASSAY PROCEDURE	23
SAMPLE PREPARATION PROCEDURE	23
TEST PROCEDURE	24
10. EVALUATION OF RESULTS	26
EXPECTED VALUES	26
11. PERFORMANCE CHARACTERISTICS	28
CROSS REACTIVITY	28
PRECISION AND REPRODUCIBILITY	28
SENSITIVITY	28
RECOVERY	29
LINEARITY	29
12. LIMITATIONS	29
13. REFERENCES	30
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	31

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Richard DM, Dawes MA, Mathias CW, Acheson A, Hill-Kapturczak N, Dougherty DM: L-Tryptophan: Basic metabolic functions, behavioural research and therapeutic indications. *International Journal of Tryptophan Research* 2009; 2 45-60.

Demisch K, Bauer J, Georgi K, Demisch L. Treatment of severe chronic insomnia with L-tryptophan: Results of a double blind cross-over study. *Pharmacopsychiatry.* 1987 Nov; 20 (6): 242-4.

Komer E, Bertha G, Flooh E, et al. Sleep inducing effect of L-tryptophane. *Er Neurol.* 1986 ; 25 (Supl 2): 107-13.

Gendall KA, Joyce PR. Meal induced changes in tryptophan:LNAAs ratio: effects on craving and binge eating. *Eat Beh* 2000; 1 (1): 53-62.

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; 8 (3): 289-95.

Breum L, Rasmussen MH, Hilsted J, Femstrom JD. Twenty-four-hour plasma tryptophan concentrations and ratio are below normal in obese subjects and are not normalised by substantial weight reduction. *Am J Clin Nutr.* 2003; 77 (5): 1112-8.

Porter RJ, Mulder RT, Joyce PR, Miller AL, Kennedy M. Tryptophan hydrolase gene (TPH1) and peripheral tryptophan level in depression. *J Affect Disord.* 1;109(1-2):209-12.

Almeida-Montes LG, Valles-Sanches V, Moreno-Aguilar J et al. Relation of serum cholesterol, lipid, serotonin and tryptophan levels to severity of depression and suicide attempts. *J Psychiatry Neurosci.* 2000 Sep; 25 (4): 371-7.

Fitzgerald P et al: tryptophan catabolism in females with irritable bowel syndrome: relationship to interferon gamma, severity of symptoms and psychiatric co-morbidity; *Neurogastroenterol Motil* 2008; 20: 1291-7

Hideki Miura et al: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress*, Volume 11, Issue 3 2008, 198-209.

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Manual

Tryptophan ELISA Kit

*For the determination of Tryptophan in human EDTA-plasma,
serum and urine*

For research use only



Valid from 17.03.2010

K 7730

96

1. INTENDED USE

The Tryptophan ELISA Kit is intended for the quantitative determination of Tryptophan in human EDTA-plasma, serum and urine. It is for research use only.

2. INTRODUCTION

Tryptophan is an essential amino acid, which can not be synthesised by humans and, therefore, must be obtained from diet. Tryptophan plays an important role in the synthesis of serotonin, as it is the only source for synthesizing serotonin (see figure 1).

Serotonin itself is released from the appendices of a widespread transmitter system with global modulating effects on different areas of life. The activity of this system depends not least on the actual tryptophan level. In some disease patterns a tryptophan deficiency is found, which is induced by an activation of the enzyme indoleamine 2,3-dioxygenase (IDO). The activation of IDO causes shifting the metabolism of tryptophan to the kynurenine pathway and its degradation (see figure 1).


Figure 1: Synthesis of serotonin and melatonin via tryptophan (left) and shifting the metabolism of tryptophan to the kynurenine-pathway via the enzyme indoleamine 2,3-dioxygenase (right)

As a consequence of that, the tryptophan level in plasma is reduced; it is not available for serotonin synthesis and in the following, for the melatonin production. As a result, fluctuations of mood, aggression, sleep and eating disorder or depression can occur.

In recent years, an increasing number of disease patterns has been detected, in which either the tryptophan level is significantly lowered due to the activation of IDO or for other reasons, or which show an adverse ratio of tryptophan-to-kynurenine levels. In the former case, the tryptophan deficiency is associated with an increase in inflammation markers. Some indications are mentioned in the following:

- In case of **irritable bowel syndrome** (IBS) interferone-gamma activates the enzyme IDO, causing a tryptophan deficiency,
- In **obese patients**, inflammation markers, which show a chronic activation of the immune system, are increased as well. IDO is also activated here, and the plasma tryptophan levels are low,
- **Stress** activates IDO likewise and facilitates the degradation of tryptophan *via* the kynurenine pathway.

In addition, gastrointestinal dysfunctions like **fructose intolerance** are associated with difficulties in absorbing tryptophan from the gut.

In all these disease patterns fluctuation of mood, sleep and eating disorder or depression are known.

The symptoms of these diseases can be eased by administration of tryptophan. However, such a treatment requires the knowledge of the tryptophan level, which can be determined appropriately with the now available test.

3. PRINCIPLE OF THE TEST

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization-reagent for tryptophan coupling. Afterwards, the treated samples and the polyclonal tryptophan-antiserum are incubated in wells of microplate coated with a tryptophan-derivative (tracer). During the incubation period, the target tryptophan in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The tryptophan in the sample displaces the antibodies out of the binding to the tracer. Therefore the concentration of the tracer-bound antibody is inverse proportional to the tryptophan concentration in the sample. During the second incubation step, a peroxidase-conjugated antibody is added to each microtiter well to detect

the anti-tryptophan antibodies. After washing away the unbound components, tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in the photometer at 450 nm. The intensity of the yellow color is inverse proportional to the tryptophan concentration in the sample; this means, high tryptophan concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. Tryptophan present in the patient samples is determined directly from this curve.

4. MATERIAL SUPPLIED

Catalog No	Content	Kit Components	Quantity
K7730MTP	PLATE	One holder with precoated strips	12 x 8 wells
K7730ST	STD	Standards	6 x 1 vial
K7730KO	CTRL 1 CTRL 2	Controls	2 x 1 vial
K7730WP	WASHBUF	Wash buffer concentrate (10 fold)	2 x 100 ml
K7730AK	AB	Anti-tryptophan antibody (lyophilized)	2 x 1 vial
K7730K	2.AB	POD antibody (concentrate)	120 µl
K7730CSP	2.ABDIL	Conjugate stabilizing buffer	24 ml
K7730RP	DERBUF	Reaction buffer	25 ml
K7730DR	DER	Derivatization reagent	2 x 6 mg
K7730LM	DMSO	Dimethylsulfoxide (DMSO)	2 ml
K7730SL	ASYBUF	Assay buffer concentrate (10-fold)	11 ml
K7730TMB	SUB	TMB substrate	25 ml
K7730AC	STOP	Stop solution	15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidist.)
- Precision pipettors and disposable tips to deliver 10-1000 μl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 10000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm
(reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μl** should be centrifuged before use to avoid loss of volume.
- The **Wash buffer concentrate (WASHBUF)** should be diluted with aqua bidist. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidist.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- **Standards (STD) and Controls (CTRL1, CTRL2)** are already diluted in the **reaction buffer (DERBUF)**. Store Standards and controls frozen at -20°C, thaw before use in the test, and re-freeze immediately after use. Standards and Controls can be re-frozen up to 3 times.
- **DMSO** could crystallize at 4°C. Dissolve the crystals at 20-25°C in a water bath.
- The **Assay buffer (ASYBUF)** should be diluted with aqua bidest. before use (11 ml concentrate + 100 ml aqua bidest.), mix well.

- The content of one vial of **derivatization reagent (DER) (6 mg)** must be dissolved in **0,75 ml DMSO**. Put the vial on a horizontal shaker for 5 min. After use, the rest of the reagent should be discarded. DER must be **prepared immediately before use**. The ELISA kit can be separated into two performances by the two DER vials. Please note: **DMSO attacks all plastics but not polypropylene products and laboratory glass**.
- The **anti-tryptophan antibody (AB)** must be dissolved in **8 ml of diluted assay buffer**. The ELISA kit can be separated into two performances by the two AB vials. **Diluted anti-tryptophan antibody is stable over a longer period. It can be stored at 2-8°C for 4 weeks**.
- The **POD antibody (2.AB)** must be diluted **1:200** in **conjugate stabilizing buffer (2.ABDIL)** (110 µl 2.AB + 22 ml 2.ABDIL). The undiluted POD antibody (2.AB) is stable at **2-8°C** until the expiry date stated on the label. **Diluted POD antibody (2.AB) is not stable over a longer period. It can be stored at 2-8°C for only 5 days**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at 2-8°C.

7. PRECAUTIONS

- For research use only.
- Human materials used in kit components were tested and found to be negative for HIV and Hepatitis B. However, for safety reasons, all kit components should be treated as if potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA-plasma, serum and urine

- Venous fasting blood and urine are suited for this test system. Blood samples are stable for one week at 2-8°C. Urine samples should be sent cooled, but they are stable for 24 h at room temperature. Otherwise they are stable for two days at 2-8°C. For longer storage, blood- and urine samples should be frozen at -20°C up to the measurement.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.

- The **EDTA-plasma-, serum- and urine** samples are analyzed without any dilution.*

*If the sample volume is less than 50 µl, a 1:1 dilution in DERBUF (reaction buffer) is recommended (25 µl sample + 25 µl DERBUF). This dilution factor should be considered for data evaluation.

Samples with visible amounts of **precipitates** should be **centrifuged** at least for 5 min at 10000 x g. The resulting supernatant is used in the assay.

- For sample preparation, a DER for coupling of tryptophan is added (details are given in the sample preparation procedure).

9. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according to the enclosed manual.

Sample preparation procedure

Samples need to be diluted in DERBUF (reaction buffer) by factor **1 + 10**. Therefore, 50 µl of sample should be diluted with 500 µl of DERBUF (reaction buffer).

Coupling of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis.

1.	Bring all reagents and samples to room temperature (18-26°C)
2.	Add 100 µl of ready to use standards (STD) , 100 µl of ready to use controls (CTRL) and 100 µl of diluted samples (SAMPLE) in the corresponding vial
3.	Add 25 µl of freshly prepared derivatization reagent (DER) into each vial (standards, controls and samples), mix well and incubate for 60 min on a shaker (180-240 rpm) at room temperature (18-26°C)

4. Afterwards add **1800 µl of diluted assay buffer** into each vial, mix well and incubate for **60 min on a shaker (180-240 rpm) at room temperature (18-26°C)**

2 x 50 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA in duplicate.

Test procedure

5. Mark the positions of **standards (STD)/controls (CTRL)/ samples (SAMPLE)** in duplicate on a protocol sheet
6. Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label
7. Wash each well **5 times** by dispensing **250 µl of diluted wash buffer** into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution
8. For the analysis in duplicate, take **2 x 50 µl of standard (STD)/control (CTRL)/samples (SAMPLE)** out of the vial and add into the respective well of the microtiter plate (PLATE)
9. Add **150 µl diluted anti-tryptophan antibody (AB)** into each well. Cover the plate tightly
10. Incubate overnight **(15-20 hours) at 2-8°C**
11. Aspirate the contents of each well. Wash each well **5 times** by dispensing **250 µl of diluted wash buffer** into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution

12.	Add 200 µl diluted POD antibody (2. AB) into each well
13.	Cover plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal shaker (180-240 rpm)
14.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution
15.	Add 200 µl of TMB substrate (SUB) into each well
16.	Incubate for 8-12 min at room temperature (18-26°C) in the dark*
17.	Add 100 µl of stop solution (STOP) into each well, mix thoroughly
18.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

*The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

10. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions, e.g. with the exact volumes for standards, controls and samples/sample treatment, standards, controls and samples are equally diluted. Therefore, **no dilution factor is required for calculation of the results.** **

**At a 1:1 dilution, the dilution factor should be considered.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Expected values

Based on internal studies of evidently healthy persons (**n=40**) a mean value of 54 µmol/l was estimated. The standard variation was 12 µmol/l.

Normal range: Serum/Plasma mean value ± 2 standard variation: 54 ± 24 µmol/l


We recommend each laboratory to develop its own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or serum pools should be analyzed with each run. Results, generated from the analysis of the control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from calibration curve. In the following an example of a calibration curve is given.

Example of calibration curve

11. PERFORMANCE CHARACTERISTICS

Cross reactivity

5-HTP (5-Hydroxy-Tryptophan)	< 0,5 %
L-Phenylalanine	< 0,1 %
L-Thyrosine	< 0,1 %

Precision and reproducibility

Intra-Assay (n=12)		
Sample	Tryptophan [μmol/l]	Standard variation (SD)
1	17,7	11,4 %
2	56,7	4,1 %

Inter-Assay (n=6)		
Sample	Tryptophan [μmol/l]	Standard variation (SD)
1	19,4	6,9 %
2	59,1	4,0 %

Sensitivity

The sensitivity was set as $B_0 + 2SD$. The zero-standard was measured 6 times.

Sample	Tryptophan mean value [OD]	Standard variation (SD)	Detection limit [μmol/l]
0	2,3	0,05 %	7

Recovery

One sample was spiked with different tryptophan concentrations and measured using this assay. The analytical recovery rate was determined by the expected and measured tryptophan levels. The expected levels were calculated as the sum of the measured tryptophan concentration in the original sample and the spiked tryptophan amount. The mean recovery rate for all concentrations was 100.5 % (n=6).

Spike [μmol/l]	Tryptophan expected [μmol/l]	Tryptophan measured [μmol/l]	Recovery [%]
0	x	x=44,8	100,0
25	44,8+25=69,8	74,7	107,0
50	44,8+50=94,8	89,5	94,4

Linearity

The linearity of the ELISA was determined by the dilution of a spiked patient sample. The mean linearity was 99.8 %.

Dilution	Measured [μmol/l]	Expected [μmol/l]	Recovery [%]
original	44,8	44,8	100,0
1+1	22,5	22,4	100,4
1+3	11,1	11,2	99,1

12. LIMITATIONS

Strong hemolytic and lipemic samples often show wrong concentrations. Do not measure hemolytic and lipemic samples.

13. REFERENCES

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydrolase-1 A218C polymorphysmus is assoziated with diagnosis, but not suicidal behavior, in borderline personality disorder. *AM J Med Genet B Neuropsychiatr Genet.* 2009 Mar 5; 150B (2): 202-8.

Richard DM, Dawes MA, Mathias CW, Acheson A, Hill-Kapturczak N, Dougherty DM: L-Tryptophan: Basic metabolic functions, behavioural research and therapeutic indications. *International Journal of Tryptophan Research* 2009; 2 45-60.

Demisch K, Bauer J, Georgi K, Demisch L. Treatment of severe chronic insomnia with L-tryptophan: Results of a double blind cross-over study. *Pharmacopsychiatry.* 1987 Nov; 20 (6): 242-4.

Komer E, Bertha G, Flooh E, et al. Sleep inducing effect of L-tryptophane. *Er Neurol.* 1986 ; 25 (Suppl 2): 107-13.

Gendall KA, Joyce PR. Meal induced changes in tryptophan:LNAAs ratio: effects on craving and binge eating. *Eat Beh* 2000; 1 (1): 53-62.

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; 8 (3): 289-95.

Breum L, Rasmussen MH, Hilsted J, Femstrom JD. Twenty-four-hour plasma tryptophan concentrations and ratio are below normal in obese subjects and are not normalised by substantial weight reduction. *Am J Clin Nutr.* 2003; 77 (5): 1112-8.

Porter RJ, Mulder RT, Joyce PR, Miller AL, Kennedy M. Tryptophan hydrolase gene (TPH1) and peripheral tryptophan level in depression. *J Affect Disord.* 1;109(1-2):209-12.

Almeida-Montes LG, Valles-Sanches V, Moreno-Aguilar J et al. Realtion of serum cholesterin, lipid, serotonin and tryptophan levels to severity of depression and suicide attempts. *J Psychiatry Neurosci.* 2000 Sep; 25 (4): 371-7.

Fitzgerald P et al: tryptophan catabolism in females with irritable bowel syndrome: relationship to interferon gamma, severity of symptoms and psychiatric co-morbidity; *Neurogastroenterol Motil* 2008; 20: 1291-7

Hideki Miura et al: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress*, Volume 11, Issue 3 2008, 198-209.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for research use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number