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Glutamat Kit

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Glutamate Kit



For the determination of glutamate in human EDTA plasma and
serum

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1. INTENDED USE	14
2. INTRODUCTION	14
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	17
TEST PROCEDURE	18
10. EVALUATION OF RESULTS	19
EXPECTED VALUES	19
11. PERFORMANCE CHARACTERISTICS	20
PRECISION AND REPRODUCIBILITY	20
SENSITIVITY	21
RECOVERY	21
LINEARITY	22
12. LIMITATIONS	22
13. REFERENCES	22
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Castellanos M, et al: High plasma glutamate concentrations are associated with infarct growth in acute ischemic stroke. Neurology. 2008 Dec 2; 71(23): 1862-8

Castillo J, Martínez F, Corredera E, Aldrey JM, Noya M: Amino acid transmitters in patients with headache during the acute phase of cerebrovascular ischemic disease. Stroke. 1995 Nov; 26(11): 2035-9

Baad-Hansen L, Cairns BE, Ernberg M, Svensson P: Effect of systemic monosodium glutamate (MSG) on headache and pericranial muscle sensitivity. Cephalalgia. 2009 Apr 28

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Manual

Glutamate Kit



For the determination of glutamate in human EDTA plasma and serum

For research use only

Valid from 09.06.2010

K 7731

1. INTENDED USE

This enzyme test is intended for the determination of glutamate in human EDTA plasma and serum. It is for research use only.

2. INTRODUCTION

Glutamic acid is a non-essential amino acid. The salts of glutamic acid are known as glutamates. Glutamate plays a central role in cellular metabolism. A significant process in amino acid degradation is transamination, in which the amino group of an amino acid is transferred to an α -keto acid, typically catalyzed by a transaminase. A very common α -keto acid is α -ketoglutarate, an intermediate in the citric acid cycle. Vice versa, transamination of α -ketoglutarate gives glutamate. The resulting α -keto acid products pyruvate and oxaloacetate are key components in fundamental processes such as glycolysis, gluconeogenesis and also the citric acid cycle.

Glutamate also plays an important role in the body's disposal of excess or waste nitrogen. Glutamate undergoes deamination, an oxidative reaction catalyzed by glutamate dehydrogenase.

In the vertebrate nervous system Glutamate is the most abundant excitatory neurotransmitter. At chemical synapses it binds to glutamate receptors, such as the NMDA receptor. Because of its role in synaptic plasticity, glutamate is involved in cognitive functions like learning and memory in the brain.

Glutamate is rapidly removed from the extracellular space by Glutamate transporters in neuronal and glial membranes. In brain injury or disease, they can work in reverse and excess glutamate can accumulate outside cells, which leads to neuronal damage and eventual cell death. This excitotoxicity due to glutamate occurs as part of the ischemic cascade and is associated with stroke and diseases like amyotrophic lateral sclerosis, lathyrism, autism, some forms of mental retardation and Alzheimer's disease.

Moreover, in the brain Glutamate also serves as the precursor for the synthesis of the inhibitory neurotransmitter GABA in GABA-ergic neurons.

Free glutamic acid is present in a wide variety of foods and is often used as a food additive and flavour enhancer in the form of its sodium salt, monosodium glutamate.

3. PRINCIPLE OF THE TEST

This assay is a photometric test intended for the determination of L-glutamate by enzymatic dehydration, in which NAD^+ is transformed to NADH.

In this reaction L-glutamate is oxidized to α -ketoglutarate by reducing NAD^+ to NADH. This reaction can be measured at 340 nm and it is proportional to the amount of oxidized L-glutamate.

4. MATERIAL SUPPLIED

Catalog No	Content	Kit Components	Quantity
K7731MTP	PLATE	Microtiter plate	12 x 8 wells
K7731ST	STD	Standards	2 x 1 vial
K7731KO	CTRL 1 CTRL 2	Controls	2 x 1 vial
K7731AP	ASYBUF	Assay buffer	16 ml
K7731RP	REABUF	Reaction buffer	2 x 1 vial
K7731EB	ENZ B	Glutamate dehydrogenase	2 x 1 vial

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Double distilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000 μl
- Foil to cover the microtiter plate
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 10000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 340 nm
- Incubator 37°C

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- **Standards (STD) and controls (CTRL1, CTRL2)** are lyophilized and must be reconstituted in **300 µl bidest. water**. After use standards and controls can be stored at -20°C and can be re-frozen up to 2 times. Re-freeze immediately after use.
- The content of one vial of **reaction buffer (REABUF)** must be dissolved in **3 ml bidest. water**, mix well. Discard any remaining quantity after use. By providing two REABUF vials the kit can be separated into two performances.
- To the content of one vial of **glutamate dehydrogenase (ENZ B) 2,6 ml assay buffer (ASYBUF)** must be added. Mix well. Discard any remaining quantity after use. By providing two ENZ B vials the kit can be separated into two performances.
- All other test reagents are stable until date of expiry (see label) when stored at 2-8°C.

7. PRECAUTIONS

- For research use only.
- Reagents should not be used beyond the expiry date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA plasma and serum

- Venous fasting blood is suited for this test system. As stability of fresh blood samples is very limited, they should be stored at 2-8°C and measured within 48 h. For longer storage samples should be stored frozen at -20°C.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- If less than 50 µl of sample is provided we recommend diluting the sample 1:2 in aqua dest. (25 µl sample + 25 µl aqua dest.). This dilution factor must be considered in data evaluation.
- **Samples** with visible amounts of **precipitates** should be **centrifuged** at least for 5 min at 10000 x g. The resulting supernatant is used in the assay.

9. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according to the enclosed manual.

Test procedure

1.	Mark the positions of standards (STD) /controls (CTRL) / samples (SAMPLE) in duplicate on a protocol sheet
2.	Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at room temperature.
3.	Add 2 x 50 µl of standards (STD) / controls (CTRL) / samples (SAMPLE) into the respective well of the microtiter plate (PLATE)
4.	Add 50 µl of bidest. water into each well.
5.	Add 100 µl assay buffer (ASYBUF) into each well
6.	Add 50 µl reaction buffer (REABUF) into each well and determine absorption immediately with an ELISA reader at 340 nm (OD _{BLANK}).
7.	Add 50 µl diluted glutamate dehydrogenase (ENZ B) into each well. Cover plate tightly.
8.	Incubate at 37°C for 15 minutes .
9.	Determine absorption at 340 nm (OD _{SAMPLE}).
10.	For analysis of obtained data see chapter 10 “evaluation of results”.

10. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, and samples, and with correct sample treatment), standards, controls, and samples are equally diluted. Therefore, **no dilution factor is required for calculation of the results.** **

**If samples are diluted 1:2, the results must be multiplied by 2.

For calculation of results subtract OD values of the blank (OD_{BLANK}) from OD values after the addition of enzyme (OD_{SAMPLE}).

To generate a standard curve, ΔOD of the standards are plotted against the standard concentrations (see quality control protocol). With the obtained slope and y-intercept glutamate concentrations of the samples can be calculated:

$$\text{glutamate concentration} = \frac{(OD_{\text{Sample}} - OD_{\text{Blank}}) - \text{intercept}}{\text{slope}} \text{ } [\mu\text{mol/l}]$$

Expected values

Based on internal studies with serum samples of evidently healthy persons (**n=24**) a mean value of 144 $\mu\text{mol/l}$ was estimated. The standard variation was 34 $\mu\text{mol/l}$.

Normal range: mean value $\pm$ 2 $\times$ standard variation: 144 $\pm$ 68 $\mu\text{mol/l}$


We recommend each laboratory to develop its own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from the calibration curve. In the following an example of a calibration curve is given.

Example of calibration curve

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n=6)		
Sample	Glutamate [µmol/l]	Standard variation (SD) [%]
1	109,6	2,1
2	223,4	1,7

Inter-Assay (n=6)		
Sample	Glutamate [μmol/l]	Standard variation (SD) [%]
1	98,2	1,3
2	15,6	1,3

Sensitivity

The sensitivity was set as $B_0 + 2SD$. The zero-standard was measured 20 times.

Sample	Glutamate mean value [OD]	2 x Standard variation (2 x SD) [%]	Detection limit [μmol/l]
0	0,02	12,6	5,9

Recovery

One sample was spiked with different glutamate concentrations and measured in this assay. The analytical recovery rate was determined by the expected and measured glutamate levels. The expected levels were calculated as the sum of the measured glutamate concentration in the original sample and the spiked glutamate amount. The mean recovery rate for all concentrations was 100,3 % (n=6).

Spike [μmol/l]	Glutamate expected [μmol/l]	Glutamate measured [μmol/l]	Recovery [%]
0	x	x=100,7	100,0
25	100,7+50=150,7	152,6	101,3
50	100,7+100=200,7	200,1	99,7

Linearity

Linearity of the test was determined by diluting a spiked patient sample. The mean linearity was 97 %.

Dilution	Measured [μmol/l]	Expected [μmol/l]	Recovery [%]
original	100,7	100,7	100,0
1+1	49,6	50,4	98,4
1+3	23,3	25,2	92,5

12. LIMITATIONS

Hemolytic and lipemic samples may give erroneous results. Do not measure hemolytic and lipemic samples.

13. REFERENCES

Ferrari A, Spaccapelo L, Pinetti D, Tacchi R, Bertolini A: Effective prophylactic treatments of migraine lower plasma glutamate levels. *Cephalalgia*. 2009 Apr; 29(4):423-9

Peres MF, Zukerman E, Senne Soares CA, Alonso EO, Santos BF, Faulhaber MH: Cerebrospinal fluid glutamate levels in chronic migraine. *Cephalalgia*. 2004 Sep; 24(9): 735-9

Castellanos M, et al: High plasma glutamate concentrations are associated with infarct growth in acute ischemic stroke. *Neurology*. 2008 Dec 2; 71(23): 1862-8

Castillo J, Martínez F, Corredera E, Aldrey JM, Noya M: Amino acid transmitters in patients with headache during the acute phase of cerebrovascular ischemic disease. *Stroke*. 1995 Nov; 26(11): 2035-9

Baad-Hansen L, Cairns BE, Ernberg M, Svensson P: Effect of systemic monosodium glutamate (MSG) on headache and pericranial muscle sensitivity. *Cephalalgia*. 2009 Apr 28

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for research use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number