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L-Arginin ELISA Kit

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1. INTENDED USE	18
2. PRINCIPLE OF THE TEST	18
3. MATERIAL SUPPLIED	19
4. MATERIAL REQUIRED BUT NOT SUPPLIED	20
5. PREPARATION AND STORAGE OF REAGENTS	20
6. PRECAUTIONS	21
7. SPECIMEN COLLECTION AND PREPARATION	22
8. ASSAY PROCEDURE	22
<i>SAMPLE PREPARATION PROCEDURE</i>	23
<i>TEST PROCEDURE</i>	23
9. EVALUATION OF RESULTS	25
<i>EXPECTED RESULTS</i>	26
10. PERFORMANCE CHARACTERISTICS	27
<i>CROSS REACTIVITY</i>	27
<i>PRECISION AND REPRODUCIBILITY</i>	27
<i>SENSITIVITY</i>	27
<i>RECOVERY</i>	28
<i>LINEARITY</i>	28
<i>CORRELATION WITH HPLC-MS</i>	29
11. LIMITATIONS	29
12. REFERENCES	29
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	30

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12. LITERATUR

- Visser M et al. (2010). The role of asymmetric dimethylarginine and arginine in the failing heart and its vasculature. Eur J Heart Fail. Oct 5.
- Neri I et al. (2010). L-Arginine supplementation in women with chronic hypertension: impact on blood pressure and maternal and neonatal complications. J Matern Fetal Neonatal Med. Oct 20.
- Saleh AI et al. (2010). Protective Effect of L-Arginine in Experimentally Induced Myocardial Ischemia: Comparison With Aspirin. J Cardiovasc Pharmacol Ther. Oct 11.
- Tripathi P et al. (2009). Oral administration of L-arginine in patients with angina or following myocardial infarction may be protective by increasing plasma superoxide dismutase and total thiols with reduction in serum cholesterol and xanthine oxidase. Oxid Med Cell Longev. Sep-Oct;2(4): 231-7.
- Bailey SJ et al. (2010). Acute L-arginine supplementation reduces the O₂ cost of moderate-intensity exercise and enhances highintensity exercise tolerance. J Appl Physiol. Aug 19.

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Manual

L-arginine ELISA Kit



For the determination of L-arginine in human EDTA plasma

For research use only

Valid from 09.12.2010

K 7733

1. INTENDED USE

The L-arginine ELISA Kit is intended for the quantitative determination of L-arginine in human EDTA-plasma. It is for research use only.

2. PRINCIPLE OF THE TEST

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for L-arginine derivatization. Afterwards, the treated samples and the polyclonal L-arginine antiserum are incubated in wells of a microtiter plate coated with a L-arginine-derivative (tracer). During the incubation period, the target L-arginine in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The L-arginine in the sample displaces the antibodies out of the binding to the tracer. Therefore, the concentration of the tracer-bound antibody is inverse proportional to the L-arginine concentration in the sample. During the second incubation step, a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-L-arginine antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the L-arginine concentration in the sample; this means, high L-arginine concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. L-arginine present in the patient samples is determined directly from this curve.

3. MATERIAL SUPPLIED

Catalog No	Content	Kit Components	Quantity
K7733MTP	PLATE	One holder with precoated strips	12 x 8 wells
K7733ST	STD	Standards (0, 12.5, 30, 75, 150, 300 μ M), diluted in sample buffer (ready to use)	6 x 1 vial
K7733KO1 K7733KO2	CTRL 1 CTRL 2	Controls, diluted in sample buffer (ready to use)	2 x 1 vial
K7733WP	WASHBUF	Wash buffer concentrate (10-fold)	2 x 100 ml
K7733PV	SAMPLEBUF	Sample buffer	45 ml
K7733AP	ASYBUF	Assay buffer concentrate (10-fold)	2 x 5 ml
K7733AK	AB	Anti-L-arginine antibody (lyophilized)	2 x 1 vial
K7733VR	ABBUF	Antibody dilution buffer concentrate (10-fold)	2 x 1 ml
K7733K	2.AB	POD antibody (concentrate)	120 μ l
K7733CSP	2.ABDIL	Conjugate stabilizing buffer	24 ml
K7733DR	DER	Derivatization reagent	2 x 12.5 mg
K7733LM	DMF	Dimethylformamide (DMF)	2 ml
K7733TMB	SUB	TMB substrate	25 ml
K7733AC	STOP	Stop solution	15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Double distilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 10000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm
(reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with aqua bidest. **1:10** before use (**100 ml WASHBUF + 900 ml aqua bidest.**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The WASHBUF is stable at **2-8°C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month.**
- **Standards (STD) and controls (CTRL1, CTRL2) are already diluted** in sample buffer (SAMPLEBUF). Store standards and controls frozen at -20°C, thaw before use in the test, and re-freeze immediately after use. Standards and controls can be re-frozen up to 3 times.

- Dissolve the content of one vial of **derivatization reagent (DER) (12.5 mg) in 750 µl DMF**. CAUTION: DMF is toxic (see chapter 6 - precautions). Put the vial on a horizontal shaker for 5 min. Dispose of any rest of the reagent after use. DER must be **prepared immediately before use**. The ELISA kit can be separated into two performances by providing two DER vials. Please note: **DMF attacks all plastics but not polypropylene products and laboratory glass**.
- Dilute the **assay buffer concentrate (ASYBUF)** with aqua bidest. **1:10** before use (**5 ml concentrate + 45 ml aqua bidest.**), mix well. The ELISA kit can be separated into two performances by providing 2 x 5 ml assay buffer. Diluted assay buffer is stable over a longer period. It can be stored at **2-8°C for 4 weeks**.
- Dilute the **antibody dilution buffer concentrate (ABBUF)** with aqua bidest. **1:10** before use (**1 ml concentrate + 9 ml aqua bidest.**), mix well. The ELISA kit can be separated into two performances by providing 2 x 1 ml antibody dilution buffer concentrate. The diluted buffer is stable over a longer period. It can be stored at **2-8°C for 4 weeks**.
- Dissolve the **anti-L-arginine antibody (AB)** in **9 ml of diluted antibody dilution buffer (ABBUF)**. The ELISA kit can be separated into two performances by providing two AB vials. Diluted anti-L-arginine antibody is stable over a longer period. It can be stored at **2-8°C for 4 weeks**.
- Dilute the **POD antibody (2.AB) 1:200** with conjugate stabilizing buffer (2.ABDIL) (**e.g. 110 µl 2.AB + 22 ml 2.ABDIL, prepare only the required amount**). The undiluted POD antibody (2.AB) is stable at **2-8°C** until the expiry date stated on the label. Diluted POD antibody (2.AB) is not stable over a longer period. It can be stored at **2-8°C for only 5 days**.
- All other test reagents are ready for use. Test reagents are stable until the expiry date (see label of test package) when stored at 2-8°C.

6. PRECAUTIONS

- For research use only.
- Human materials used in kit components were tested and found to be negative for HIV and Hepatitis B. However, for safety reasons all kit components should be treated as if potentially infectious.

- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- DMF is toxic. It may cause harm to the unborn child and is harmful by inhalation and in contact with skin. Work under hood and apply preventive skin protection. In case of skin or eye contact flush with plenty of water and get medical attention.
- Reagents should not be used beyond the expiration date shown on kit label.

7. SPECIMEN COLLECTION AND PREPARATION

EDTA plasma

- Venous fasting blood is suited for this test system. Blood samples are stable for one week at 2-8°C. For longer storage, blood samples should be frozen at -20°C.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- The EDTA plasma samples are diluted for derivatization.

Samples with visible amounts of **precipitates** should be **centrifuged** at least for 5 min at 10000 x g. The resulting supernatant is used in the assay.

- For sample preparation, a DER (derivatization reagent) for derivatization of L-arginine is added (details are given in the sample preparation procedure).

8. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according to the enclosed manual.

Sample preparation procedure

Dilute **EDTA plasma samples** with reaction buffer by **factor 1:40**, i.e. **25 µl sample + 975 µl sample buffer (SAMPLEBUF)**.

Derivatization of standards (STD), controls (CTRL) and diluted samples (SAMPLE) is carried out in single analysis.

1.	Bring all reagents and samples to room temperature (18-26°C).
2.	Add 100 µl of ready to use standards (STD) , 100 µl of ready to use controls (CTRL) and 100 µl of diluted samples (SAMPLE) in the corresponding vial.
3.	Add 25 µl of freshly prepared derivatization reagent (DER) into each vial (standards, controls and samples), mix well and incubate for 45 min on a shaker (180-240 rpm) at room temperature (18-26°C) .
4.	Afterwards add 1250 µl of diluted assay buffer (ASYBUF) into each vial, mix well and incubate for 45 min on a shaker (180-240 rpm) at room temperature (18-26°C) .

2 x 50 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5.	Mark the positions of standards (STD)/ controls (CTRL)/ samples (SAMPLE) in duplicate on a protocol sheet.
6.	Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
7.	Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate (PLATE) should be firmly tapped on absorbent paper to remove excess solution.

8.	For the analysis in duplicate, take 2 x 50 µl of standards (STD) / controls (CTRL) / samples (SAMPLE) out of the vial and add into the respective wells of the microtiter plate (PLATE).
9.	Add 150 µl diluted anti-L-arginine antibody (AB) into each well. Cover the plate tightly.
10.	Incubate overnight (15-20 hours) at 2-8°C .
11.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.
12.	Add 200 µl diluted POD antibody (2. AB) into each well.
13.	Cover plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal shaker (180-240 rpm).
14.	Aspirate the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.
15.	Add 200 µl of TMB substrate (SUB) into each well.
16.	Incubate for 12-16 min at room temperature (18-26°C) in the dark*.
17.	Add 100 µl of stop solution (STOP) into each well, mix thoroughly.
18.	Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e.g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

9. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, samples, and with correct sample treatment), standards, controls and samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results.**

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

1. *4-parameter-algorithm*

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).

2. *Point-to-point-calculation*

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. *Spline-algorithm*


We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).

Plausibility of the measured pairs of values should be examined before automatically evaluating the results. If this option is not available within the used program, the pairs of values should be controlled manually.

Controls

Control samples or plasma pools should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from the calibration curve. In the following, an example of a calibration curve is given. Do not use it for the calculation of your results.

Example of calibration curve***Expected results***

Based on internal studies with plasma samples of evidently healthy persons (n=40) a mean value of 84.3 $\mu\text{mol/l}$ was calculated. The standard deviation was 25.3 $\mu\text{mol/l}$.

Plasma mean value \pm 2 x standard deviation: **84.3 \pm 50.6 $\mu\text{mol/l}$**

Normal range: **33.7 – 134.9 $\mu\text{mol/l}$**

We recommend each laboratory to develop its own normal range. The values mentioned above are only for orientation and can deviate from other published data.

10. PERFORMANCE CHARACTERISTICS

Cross reactivity

ADMA	< 0.06%
SDMA	< 0.006%

Precision and reproducibility

EDTA plasma:

Intra-Assay (n=8)		
sample	L-arginine [μmol/l]	coefficient of variation (CV) [%]
1	51.4	8.5
2	94.4	8.9

Inter-Assay (n=6)		
sample	L-arginine [μmol/l]	coefficient of variation (CV) [%]
1	46.4	8.6
2	95.4	3.6

Sensitivity

The sensitivity was set as $B_0 + 2SD$. The zero-standard was measured 32 times.

sample	L-arginine mean value [OD]	2 x standard deviation (2 x SD)	detection limit [μmol/l]
zero-standard	1.7	0.1	3.0

Recovery

One sample was spiked with different L-arginine concentrations and measured using this assay. The analytical recovery rate was determined by the expected and measured L-arginine levels. The expected levels were calculated as the sum of the measured L-arginine concentration in the original sample and the spiked L-arginine amount. The mean recovery rate for all concentrations was 102.2 % (n=10).

EDTA plasma:

spike [μmol/l]	L-arginine measured [μmol/l]	L-arginine expected [μmol/l]	recovery [%]
0	x = 49.5	x	100.0
50	97.0	49.5+50 = 99.5	97.5
100	162.9	49.5+100 = 149.5	109.0

Linearity

The linearity of the ELISA was determined by the dilution of a spiked patient sample. The mean linearity was 100.5 % (n=6).


EDTA-plasma

dilution	measured [μmol/l]	expected [μmol/l]	recovery [%]
original	97.0	97.0	100.0
3:4	66.7	72.8	91.6
2:3	64.6	64.7	99.8
1:2	49.5	48.5	102.1
1:3	35.2	32.3	109.0

Correlation with HPLC-MS

14 samples were measured with this ELISA and HPLC-MS. The correlation was $r = 0.982$.

HPLC-MS vs. ELISA

11. LIMITATIONS

Hemolytic and lipemic samples may give erroneous results. Do not measure hemolytic and lipemic samples.

12. REFERENCES

- Visser M et al. (2010). The role of asymmetric dimethylarginine and arginine in the failing heart and its vasculature. Eur J Heart Fail. Oct 5.
- Neri I et al. (2010). L-Arginine supplementation in women with chronic hypertension: impact on blood pressure and maternal and neonatal complications. J Matern Fetal Neonatal Med. Oct 20.
- Saleh Al et al. (2010). Protective Effect of L-Arginine in Experimentally Induced Myocardial Ischemia: Comparison With Aspirin. J Cardiovasc Pharmacol Ther. Oct 11.

- Tripathi P et al. (2009). Oral administration of L-arginine in patients with angina or following myocardial infarction may be protective by increasing plasma superoxide dismutase and total thiols with reduction in serum cholesterol and xanthine oxidase. *Oxid Med Cell Longev*. Sep-Oct;2(4): 231-7.
- Bailey SJ et al. (2010). Acute L-arginine supplementation reduces the O₂ cost of moderate-intensity exercise and enhances highintensity exercise tolerance. *J Appl Physiol*. Aug 19.

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for research use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number