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Nitrotyrosin ELISA

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
Nitrotyrosine ELISA

For the in vitro determination of Nitrotyrosine in human sample material, e.g. stool

Valid from 23.10.2009

K 7824

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of Nitrotyrosine in human sample material, e.g. stool. It is for *in vitro* use only.

2. INTRODUCTION

Nitrotyrosine is the nitrated form of the amino acid tyrosine. The accumulation of protein bound nitrotyrosine is associated with cardiovascular diseases that are based on inflammatory processes (e.g., atherosclerosis, myocardial infarction, diabetic vasculopathy, hypertension, or coronary heart diseases). A growing number of studies have also associated the accumulation of nitrotyrosine with neurological diseases (Alzheimer´s disease, Parkinson´s disease, multiple sclerosis, stroke). With treatment of some of the associated diseases the levels of nitrated tyrosines have been shown to decrease, so nitrotyrosine has been stated to be a marker of nitrosative stress.

During inflammatory processes, large amounts of nitric oxide ($\bullet\text{NO}$) are locally released from L-arginine. This reaction is catalyzed by the enzyme NO-synthase (NOS). Other causes for the increased $\bullet\text{NO}$ production are exposure to chemicals or heavy metals, drugs, nicotine, or physical and psychological stress, as well as extraordinary physical strain with increased oxygen consumption.

In high concentrations, $\bullet\text{NO}$ that is not trapped by mitochondrial superoxide dismutase (SOD) reacts with superoxide ($\text{O}_2\bullet^-$) to form peroxynitrite (ONOO^-). Peroxynitrite is implicated as a key oxidant species in several pathologies and is known to be cytotoxic (nitrosative stress).

Peroxynitrite is highly reactive and shows a high affinity to aromatic amino acids, e.g., to the phenolic ring of tyrosine. The nitration of tyrosine in general is a natural process within the post-translational protein modification.

Nitrotyrosine is a stable product and might be seen as a correlate of peroxynitrite production, and its accumulation in cells and tissues is a marker of oxidative stress and nitrosative stress, respectively (Ischiropoulos 2008).

Indications

- Cardiovascular diseases
- Neurological diseases
- Thyroid disturbances
- Blockade of biochemical pathways
- Mitochondriopathy

Consequences of nitrosative stress

- Modification of lipids and proteins (for example structural proteins in mitochondria)
- Inhibition of respiratory chain enzymes in the mitochondria
- Glutamate overload
- Disturbances in ion channels
- Calcium overload
- Initiation of apoptosis processes

3. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 7824MTP	PLATE	Microtiter plate, precoated	12 x 8 wells
K 7824WP	WASHBUF	ELISA wash concentrate 10x	1 x 100 ml
K 7824K	CONJ	HRP-antibody (concentrate)	1 x 150 µl
K 7824ST	STDKONZ	Standard (concentrate)	1 x 240 µl
K 7824AP	ASYBUF	Assay buffer	1x 50 ml
K 7824KO1	CTRL POS	Positive control	4 x 250 µl
K 7824KO2	CTRL NEG	Negative control	4 x 250 µl
K 7824TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 7824AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Standard laboratory polypropylene reaction vessels (1.5 ml)
- Standard laboratory reaction vessel (15 ml)
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Vortex mixer
- Bidistilled water (aqua bidest.)
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- Reagents with a volume less than 100 μL should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C using a water bath before dilution. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The controls (**CTRL POS** and **CTRL NEG**) must be reconstituted with 500 μl of **ASYBUF** and used in the assay without any further dilution.
- The **CONJ** (HRP-antibody) must be diluted **1:100** with diluted wash buffer. It should be freshly prepared for each run. The remainder of the diluted conjugate should be discarded.
- For Standard 1, pipet **50 μl** of **STD** (Standard stock solution) into 450 μl **ASYBUF**. It corresponds to a concentration of 1250 nM nitrotyrosine.

A standard curve is prepared in **1 : 2.5** dilution steps according to the following pipetting scheme.

Standard 1 : 1250 nM	50 μl STD (Standard stock) + 450 μl ASYBUF
Standard 2 : 500 nM	200 μl Standard 1 + 300 μl ASYBUF
Standard 3: 200 nM	200 μl Standard 2 + 300 μl ASYBUF
Standard 4: 80 nM	200 μl Standard 3 + 300 μl ASYBUF
Standard 5: 32 nM	200 μl Standard 4 + 300 μl ASYBUF
Standard 6: 0 nM	ASYBUF

Important: Polypropylene reaction vessels should be used for the standard dilution series.

The standard curve should be freshly prepared for each run. The remainder should be discarded.

- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

6. SAMPLE PREPARATION

Liquid samples are diluted with ASYBUF (assay buffer). A working dilution of 1:5 is advised but optimal dilution should be determined by experimentation.

Serum, Plasma: 1:5 in ASYBUF (assay buffer)

Extraction of stool samples

The assay buffer is used as a sample extraction buffer. We recommend the following sample preparation:

1. Use of the stool sample preparation kit for dosing a 100 mg stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). The stool sample must be suspended in 5 ml ASYBUF (assay buffer). This corresponds to a dilution factor of 1:50.
2. Afterwards, mix stool sample and buffer; vortex for at least 30 sec. depending on the stool consistency.
3. Thereafter, transfer ca. 1 ml of the stool suspension into a closable vial, e.g. Eppendorf-tube and centrifuge for 5 minutes at 13000 rpm.

The stool suspension is not stable. It is recommended to use a fresh sample for each run.

After the centrifugation, the supernatant is further diluted with ASYBUF (assay buffer). We recommend a dilution of 1:2 to 1:10, which corresponds to a **final dilution of 1:100 – 1:500** of the stool sample.

For analysis, pipette **100 µl of the final dilution** per well.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the “sandwich” technique with two polyclonal antibodies against nitrated proteins.

Standards, controls and diluted samples which are assayed for nitrotyrosine are added into the wells of a micro plate coated with polyclonal rabbit anti-nitrotyrosine antibody. During the first incubation step, nitrated proteins are bound by the immobilized primary antibody. Then a peroxidase-conjugated polyclonal rabbit anti- nitrotyrosine antibody is added into each microtiter well and a “sandwich” of

primary antibody - nitrated protein – peroxidase-conjugate

is formed. Tetramethylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of nitrotyrosine. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. standard concentration is generated, using the values obtained from the standard.

Test procedure

1. Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well
2. Mark the positions of STD (Standard), CTRL POS (positive control), CTRL NEG (negative control) and SAMPLE (sample) on a protocol sheet
3. Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label
4. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper

5.	Add 2 x 100 µL of the prepared STD (Standard), CTRL POS (positive control), CTRL NEG (negative control) and SAMPLE (sample) in duplicate into respective well
6.	Cover plate or strips with foil tightly and incubate for 2.5 h at room temperature (18 - 26°C) on the horizontal shaker
7.	Discard the contents of each well. Wash 5 times by dispensing 250 µL of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
8.	Pipette 100 µL of diluted CONJ (HRP-antibody) into each well
9.	Cover plate or strips with foil tightly and incubate for 1h at room temperature (18 - 26°C) on the horizontal shaker.
10.	Discard the contents of each well. Wash 5 times by dispensing 250 µL of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
11.	Add 100 µL of SUB (TMB substrate) into each well
12.	Incubate for 10-20 min at room temperature (18-26°C) in the dark.
13.	Add 100 µL of STOP (stop solution) into each well, mix thoroughly in a microtiter plate reader
14.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

*The intensity of the colour change is temperature sensitive. We recommend to observe the colour change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-Parameter-algorithm".

1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline-algorithm

We recommend a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Samples

The obtained nitrotyrosine concentration must be multiplied by the corresponding dilution factor.

Controls

The nitrotyrosine concentration can be read directly from the calibration curve. The concentration range is given on the enclosed data sheet specification.

9. LIMITATIONS

Samples with nitrotyrosine concentrations **outside the standard curve range** should be diluted further with ASYBUF (assay buffer) to obtain readings within the standard curve and re-assayed.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

It is recommended that each laboratory should determine its own normal range.

11. PERFORMANCE CHARACTERISTICS

Sensitivity

Detection limit: 34,1 nM

The sensitivity limit was set as $B_0 + 3 \text{ SD}$. The Zero-standard was measured 8 times.

Cross reactivity

Albumin treated with nitrating agents: 100 %

Albumin non-treated: < 0,1 %

12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulphuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

15. REFERENCES

Gonsette RE (2008) Neurodegeneration in multiple sclerosis : The role of oxidative stress and excitotoxicity. J Neurol Sci, Vol 274, Issue 1-2, 48-53

Ischiropoulos H (2008) Protein tyrosine nitration – An update. Arch Biochem Biophys Oct 30

Peluffe G, Radi R (2007) Biochemistry of protein tyrosine nitration in cardiovascular pathology. Cardiovasc Res. Jul 15 : 75(2) :291-302

Souza JM et al. (2008) Protein tyrosine nitration- functional alteration or just a biomarker ? Free Radic Biol Med. Aug 15 ; 45 (4) :357-356

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		For research use only