# Biotin ELISA Kit

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<b>INHALTSVERZEICHNIS/TABLE OF CONTENT</b>	<b>SEITE/PAGE</b>
<b>1. INTENDED USE</b>	<b>19</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>19</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>19</b>
<b>4. MATERIAL SUPPLIED</b>	<b>20</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>21</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>21</b>
<b>7. PRECAUTIONS</b>	<b>22</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>22</b>
<b>9. ASSAY PROCEDURE</b>	<b>23</b>
PROCEDURAL NOTES	23
TEST PROCEDURE	24
<b>10. RESULTS</b>	<b>25</b>
<b>11. LIMITATIONS</b>	<b>26</b>
<b>12. QUALITY CONTROL</b>	<b>26</b>
EXPECTED VALUES	26
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>27</b>
SPECIFICITY	27
SENSITIVITY	27
RECOVERY	27
INTRA- AND INTERASSAY VARIATION	28
VARIATION IN CALIBRATION CURVE	28
DILUTION LINEARITY	29
<b>14. REFERENCES</b>	<b>30</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>31</b>

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17.10.2006 10012002.doc


# Biotin ELISA Kit

*For the in vitro determination of Biotin (Vitamin H) in Serum,  
Plasma, Drugs and Food*

Gültig ab/valid from 17.10.2006

K 8140

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Biotin** in serum, plasma and foods. For research use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

In the late 1950s Lynen (1) and Wakil (2) reported the discovery of covalently bound biotin with a coenzyme function. Biotin is a growth factor present in minute amounts in every living cell. It plays an indispensable role in numerous naturally occurring carboxylation reactions. The daily biotin requirements for adults lie between 100 and 200 mg. Renewed interest in the role of biotin in human nutrition and therapy stems from the accumulating reports of various biotin-responsive syndromes (3-5). A variety of methods for the assay of biotin are currently available (6-11). However most of them are either cumbersome to perform or based on the use of radioactive materials. This assay is a convenient reciprocal assay for biotin. It is extremely sensitive, easy to perform and doesn't use radiolabels.

## 3. PRINCIPLE OF THE TEST

The test is a competitive enzyme-linked assay for biotin determination in human serum. Samples, standards and controls are incubated with streptavidin-enzyme conjugate for 15 minutes at 18-25 °C. Then 100 µl of each incubation mixture is transferred in duplicate into microtiter plate wells coated with biotin-albumin conjugate. After a second incubation step of 30 minutes at 18-25 °C, the enzyme substrate TMB is added. Finally, the reaction is terminated by an acidic stop solution causing a color change from blue to yellow. The color intensity is inversely proportional to the biotin concentration. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Biotin present in the patient samples is determined from this curve.

**4. MATERIAL SUPPLIED**

Cat. No	Content	Kit Components	Quantity
K 8140MTP	PLATE	one holder with precoated strips	96
K 8140WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 25 ml
K 8140VP	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 15 ml
K 8140KP	CONJBUF	Conjugate dilution buffer, ready-to-use	1 x 25 ml
K 8140ST	STD	Calibrators, ready-to-use (red cap) (0; 37.5; 75; 150; 300; 600 ng/l)	6 x 1 ml
K 8140K	CONJ	Conjugate, horseradish peroxidase labeled streptavidin (green cap)	1 x 0.3 ml
K 8140KO	CTRL	Controls, ready-to-use (yellow cap)	2 x 1 ml
K 8140TMB	SUB	TMB Substrate (Tetramethylbenzidin), ready-to-use	1 x 12 ml
K 6750AC	STOP	ELISA Stop solution, ready-to-use	1 x 12 ml

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Distilled or deionized water.
- Precision pipettes with disposable tips to deliver in the ranges 2-20 µl,
- 20-200 µl and 200-1000 µl.
- Clean disposable test tubes.
- Range of standard, clean volumetric laboratory glassware.
- Absorbent paper towels.
- Vortex mixer.
- Microtiter plate reader with 450 nm filter.

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua dest. **1:10** before use (add 450 ml aqua dest. to 50 ml concentrate). Crystals could occur due to high salt concentration. The crystals must be resuspended **before dilution of the buffer solutions** using a water bath at 37°C. The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions can be stored at 2-8°C for 2 weeks.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **Conjugate** (CONJ) must be diluted **1:100** in conjugate dilution buffer (CONJBUF). The antibody is stable at 18-25 °C in the dark for 8 hours.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2 – 8 °C**.

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution (STOP) consists of diluted HCl. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### Serum and plasma

Human serum can directly be processed but it should be free from insoluble particles. Remove insoluble particles by centrifugation. Serum samples can be stored at 2-8 °C for one day. For long time storage serum samples should be frozen. Repeated freezing and thawing should be avoided. Serum samples with expected high biotin concentrations (exceeding 600 ng/liter) should be properly diluted with the provided sample diluent (SAMPLEBUF).

### Milk

Milk can be assayed without pre-treatment but should be diluted 1:100 in sample dilution buffer in the following way: first, add 50 µl milk to 450 µl sample dilution buffer (SAMPLEBUF) in a test tube, vortex and then add 50 µl of this solution to 450 µl sample dilution buffer in a second test tube.

***Instructions for washing***

Wash each well 5 times by dispensing 0.3 ml of diluted wash solution into each well. Aspirate completely the contents of each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.

Alternatively, automatic washers may be used. Ensure that automatic washers are functioning correctly, i.e. that the contents of the wells are aspirated completely and the bottoms of the wells are not scratched.

**9. ASSAY PROCEDURE*****Procedural notes***

- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variation of the test procedure, that is not coordinated with the producer, may influence the test results. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The incubation time (Point 4) must be 30 min. The pipetting time of the samples should not take longer than 10 min.
- The assay should always be performed according the enclosed manual.

## Test procedure

1. Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label. Allow all reagents to come to room temperature (18-25 °C) prior to use.
2. Add **150 µl** of STD/SAMPLE/CTRL (Standard/Sample/Control) into individual test tubes.
3. Add **150 µl** of diluted streptavidin-enzyme conjugate solution into each test tube, vortex-mix.
4. Incubate for **30 minutes** at 18-25 °C.
5. Transfer **100 µl** of the content of each tube in duplicate into individual microtiter wells.
6. Incubate for **30 minutes** at 18-25 °C.
7. **Wash wells 5 times** (see instructions for washing).
8. Pipet **100 µl** of the TMB substrate (SUB) solution into all wells.
9. Incubate for **10-30 minutes** at 18-25 °C.\*
10. Stop the reaction by adding **100 µl** of hydrochloric acid (STOP) to each well and mix carefully.
11. Read OD on ELISA reader at **450 nm** within 30 minutes after adding the stop solution. For dual wavelength analysis, 620 nm, 655 nm or 690 nm should be used as a reference wavelength.

\*The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good color differentiation. For the *in vitro* *biotin* test it takes approximately 10 minutes.

## 10. RESULTS

Calculate the average absorbance values for each set of duplicate standards and samples as a percentage of the absorbance value of the zero standard. Duplicates should be within 20 percent of the mean.

$$\frac{\text{Absorbance standard/sample}}{\text{Absorbance zero standard}} \times 100\% = \% \text{ Absorbance}$$

Draw a standard curve by plotting the log biotin concentration on the ordinate against the % absorbance for each standard concentration on the abscissa. Determine the concentration biotin (ng/liter) for each sample. The concentration read from the standard curve must be further multiplied by the corresponding dilution factor.

A representative standard curve is shown in figure 1. This curve cannot be used to evaluate test results. Each laboratory must prepare a standard curve for each test run.

Table 1.

Biotin (ng/liter)	Absorbance (450 nm)	% Absorbance
600	0.210-0.210	10.5
300	0.346-0.347	17.4
150	0.783-0.777	39.1
75	1.574-1.566	78.6
37.5	1.876-1.880	94.0
0	2.041-1.952	100


Figure 1

## 11. LIMITATIONS

Samples with Biotin levels greater than the highest calibrator value should be diluted and re-assayed.

## 12. QUALITY CONTROL

### *Expected values*

Healthy:  $\geq 200$  ng/l

Optimal Status: 100 - 200 ng/l

Vitamin deficiency:  $< 100$  ng/l

## 13. PERFORMANCE CHARACTERISTICS

### *Specificity*

The specificity of the *BIOTIN* assay was established by analyzing the cross-reactivity with related compounds. The values represent the cross-reactivity as a percentage of the biotin reactivity; biotin-4-amidobenzoic acid - 81%, biocytin - 106%.

### *Sensitivity*

The limit of detection was defined as the lowest biotin concentration distinguishable from the zero standard and calculated as the concentration causing an absorbance value lower than the absorbance of the zero standard minus 3 SD, according to the recommendation of the International Federation of Clinical Chemistry (12). The detection limit was determined to be 15 ng biotin/liter. The effective dose for 50% of maximal response was 100 ng biotin/liter.

### *Recovery*

Human serum with 86 ng/liter endogenous biotin was spiked with variable amounts of biotin and assayed with the *in VITRO TEST BIOTIN* (table 2).

Table 2.

Spiked with (ng/liter) biotin	Biotin expected (ng/liter)	Biotin measured (ng/liter)	Recovery (%)
250	336	366	109
125	211	213	101
62.5	148.5	141	95
31.25	117.25	110	94

The overall recovery is 100% ± 7%.

### *Intra- and interassay variation*

The intra- and interassay variation was respectively 1.5% and 4.9%.

### *Variation in calibration curve*

The precision profile was generated by preparing five replicates of each of the standards and calculating the standard variation and percent coefficient of variation (CV) at each concentration (figure 2). The average CV between the standards was 1.8 %.


Figure 2

## Dilution linearity

The dilution linearity profile was generated by diluting serum, with a normal amount of endogenous biotin, with sample diluent. The assay was carried out as described.


Figure 3

## 14. REFERENCES

1. Lynen, F., Knappe, J., Lorch, E., Jutting, G. and Ringelmann, E. (1959) *Angew. Chem.* **71**, 481.
2. Wakil, S. J. and Gibson, D.M. (1960) *Biochim. Biophys. Acta* **41**, 122.
3. Bonjour, J.P. (1977) *Int. J. Vit. Nutr. Res.* **47**, 107.
4. Dakshinamurti, K. and Bhagavan, H.N. eds. (1985) Biotin, Vol. 447, pp. 1-441, *Ann. N.Y. Acad. Sci.*, New York
5. Watanabe, T. (1983) *J. Nutr.* **113**, 574.
6. Green, N.M. (1970) in *Methods in Enzymology* (McCormick, D.B. and Wright, L.D., eds), Vol. 18, pp. 383-385, Academic press, New York.
7. Hood, R.L. (1979) in *Methods in Enzymology* (McCormick, D.B. and Wright, L.D., eds), Vol. 62, pp. 279-283, Academic Press, New York.
8. Dakshinamurti, K. and Allan, L. in *Methods in Enzymology* (McCormick, D.B. and Wright, L.D., eds), Vol. 62, pp. 284-287, Academic press, New York. 11 January 2006. Green, N.M. (1970) in *Methods in Enzymology* (McCormick, D.B. and Wright, L.D., eds), Vol. 18, pp. 418-424, Academic press, New York.
9. Lin, H.J. and Kirsch, J.F. (1979) in *Methods in Enzymology* (McCormick, D.B. and Wright, L.D., eds), Vol. 62, pp. 287-289, Academic press, New York.
10. Al-Hakiem, M.H.H., Landon, J., Smith, D.S. and Nargessi, R.D. (1981) *Anal. Biochem.* **116**, 264.
11. Büttner, J., Borth, R., Boutwell, J.H., Broughton, P.M.G. and Bowyer, R.C. (1979) Quality control in clinical chemistry. *Clin. Chim. Acta* **98**, 145F-162F.
12. Livaniou, E., Evangelatos, G.P. and Ithakissios, D.S. (1987) *Clin. Chem.* **33**, 1983-1988.

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- For *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

17.10.2006\_050201.doc