anti-Gliadin IgA ELISA

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CONTENT	PAGE
1. INTENDED USE	15
2. INTRODUCTION	15
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	16
4. MATERIAL SUPPLIED	17
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	18
7. PRECAUTIONS	19
8. SPECIMEN COLLECTION AND PREPARATION	19
9. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
TEST PROCEDURE	20
10. LIMITATIONS	21
11. RESULTS	21
12. QUALITY CONTROL	21
EXPECTED VALUES	21
13. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
14. TECHNICAL HINTS	23
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
16. REFERENCES	23

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Bestellnummer

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
anti-Gliadin IgA ELISA

For the in vitro determination of IgA anti-gliadin antibodies in serum

Valid from 09.07.2010

K 9310

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the determination of **IgA anti-gliadin antibodies** in serum. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

Celiac disease (CD), also known as gluten sensitive enteropathy (GSE) or celiac sprue, is an immune-mediated disease, triggered by gluten proteins in wheat, rye and barley. It is mainly manifested as chronic intestinal malabsorption during childhood. Celiac disease is particularly associated with T cell lymphomas of the small intestine, which occur in 6 of 10 celiac sprue patients. Thus, early diagnosis of celiac disease is important, as people can avoid long-term health complications of celiac disease by following a gluten-free diet.

The determination of IgA anti-gliadin antibodies is a sensitive tool for screening, therapy control and monitoring remission of celiac disease along with distal duodenal biopsies, the gold standard for diagnosis.


Indication

- Auto immune-mediated diseases
- Food intolerance
- See our proposal for stepwise diagnosis of gluten sensitivity on page 16.

3. PRINCIPLE OF THE TEST

The wells of the microtiter plate are coated with the antigen gliadin. In a first incubation step, the anti-gliadin antibodies are bound to the coated antigen. To remove all unbound substances, a washing step is carried out.

In a further incubation step, a peroxidase-labeled rabbit IgA anti-gliadin antibody is added. After another washing step to remove all unbound substances, the solid phase is incubated with the peroxidase substrate, tetramethylbenzidine (TMB). Finally, an acidic stop solution is added to terminate the enzymatic reaction, whereby the colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the IgA anti-gliadin antibody concentration. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standards. IgA anti-gliadin antibodies present in the samples are determined directly from this curve.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 9310MTP	PLATE	One holder with pre-coated strips	8 x 12 wells
K 9310WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 9310K	CONJ	Conjugate (rabbit-anti-IgA antibody, peroxidase labeled)	1 x 200 µl
K 9310ST	STD	Gliadin Standards, lyophilized	4 x 6 vials
K 9310KO1	CTRL	Control, lyophilized	4 x 1 vial
K 9310KO2	CTRL	Control, lyophilized	4 x 1 vial
K 9310TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 9310AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a **volume less than 100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be re-dissolved at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The **CONJ** (conjugate, peroxidase-labeled) must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The **undiluted conjugate** is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable at **2-8 °C** until expiry date stated on the label. Each of the lyophilized standards and controls must be reconstituted with **500 µl of bidistilled water** (aqua bidest.). Allow the vial content to dissolve for 10 minutes at room temperature, and mix thoroughly by gentle inversion to insure complete reconstitution. **Diluted standards and controls are not stable and can not be stored.**
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum/Plasma

Serum must be diluted **1:250** in ELISA wash buffer in two steps.

For example:

Dilution step 1: 10 µl serum + 90 µl wash buffer

Dilution step 2: 40 µl dilution 1 + 960 µl wash buffer

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

Wash the pre-coated microtiter plate **5 x with 250 µl** of wash buffer before use. Carry out the tests in duplicate.

1. Add **100 µl** of **STD** (standards), **CTRL** (controls) and **diluted samples**.
2. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
3. Decant the content of the plate and wash the wells **5x with 250 µl** of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel.
4. Add **100 µl** of diluted **CONJ** (conjugate).
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5x with 250 µl** of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel.
7. Add **100 µl** of **SUB** (TMB substrate).
8. Incubate for **10 - 15 minutes** at room temperature.
9. Add **50 µl** of **STOP** (ELISA stop solution) and mix shortly.
10. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

10. LIMITATIONS

Samples with IgA anti-gliadin levels higher than the highest standard should be further diluted with wash buffer and re-assayed.

11. RESULTS

Calculation

Since the sample dilution of 1:250 is already considered in the calibration curve, the dilution factor is 1.

If a sample dilution different from 1:250 has been used, a corresponding dilution factor must be considered, e.g., at a dilution of 1000 the results must be multiplied by 4.

12. QUALITY CONTROL

Controls should be analyzed with each run. Results, generated from the analysis of controls, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality controls are outside the acceptable limits.

Expected values

Preliminary normal value

Serum > 22 AU/ml positive

Serum < 16 AU/ml negative

Results between 16 - 22 AU/ml are in the grey range.

We recommend each laboratory to establish its own normal range.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay-Variation

The intra-assay variation of the assay was calculated from 64 determinations of one sample by one person.

Intra-Assay CV n = 64

Sample	IgA anti-gliadin mean value [AU/ml]	Intra-Assay CV [%]
1	133,41	4,9

Inter-Assay-Variation

The inter-assay variation of the assay was calculated from data on 2 samples obtained in 17 determinations by different persons on different days.

Inter-Assay CV n = 17

Sample	IgA anti-gliadin mean value [AU/ml]	Inter-Assay CV [%]
1	133,53	7,8
2	5,04	18,4

14. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in-vitro* diagnostics only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

16. REFERENCES

1. Dieterich et al., 1997; Nat. Med. 7 (3), 797
2. Rieken et al., 1998; Dtsch. med. Wschr. 123, 1454
3. Green et al., 1998; Clin. Persp. Gastroenterol., November, 133
4. Mothes, Th., 1997; Münsch. Med. Wschr. 139, 111
5. Sardy et al., 1999; Clin. Chem 45, 2142

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number