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anti-Gliadin slgA ELISA

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Additional information about our products is available on our homepage

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Table of content

1. INTENDED USE	17
2. INTRODUCTION	17
3. PRINCIPLE OF THE TEST	17
4. MATERIAL SUPPLIED	18
4. MATERIAL SUPPLIED	19
5. MATERIAL REQUIRED BUT NOT SUPPLIED	19
6. PREPARATION AND STORAGE OF REAGENTS	20
7. PRECAUTIONS	20
8. SPECIMEN COLLECTION AND PREPARATION	21
9. ASSAY PROCEDURE	22
PROCEDURAL NOTES	22
TEST PROCEDURE	23
10. RESULTS	24
11. QUALITY CONTROL	24
EXPECTED VALUES	24
12. PERFORMANCE CHARACTERISTICS	25
PRECISION AND REPRODUCIBILITY	25
13. REFERENCES	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26

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Manual



anti-Gliadin slgA ELISA

For the in vitro determination of anti-Gliadin slgA antibodies in stool

Valid from 25.03.2011

K 9311

1. INTENDED USE

The Assay is intended for the quantitative determination of **anti-gliadin-slgA antibodies** in stool. For *in vitro* diagnostic use only.

2. INTRODUCTION

The celiac/coeliac disease is caused by the gliadin fraction of wheat gluten and similar proteins of rye and barley. The disease is mainly manifested as chronic digestion insufficiency in children or young adults. In addition, patients with celiac disease have a greatly increased risk of developing malignant T-cell lymphoma of the small bowel, as T-cell lymphoma was found in 6 of 10 patients with coeliac sprue. Early diagnosis of celiac disease is important, because there is evidence that a gluten-free diet might help to prevent complications and malignancies (including intestinal lymphoma).


The measurement of **anti-gliadin-slgA-antibodies** may be useful as a screening criterion before jejunal biopsy, the „Gold Standard“, and for the monitoring of the gluten-free diet treatment.

Indication

- Celiac disease
- Food intolerance
- See also our proposal for stepwise diagnosis of gluten sensitivity on page 18.

3. PRINCIPLE OF THE TEST

The antigen gliadin is immobilized on the microtiter plate. During the first incubation step, the human anti-gliadin-slgA antibodies in the samples are bound by the immobilized antigen. After a washing step, the anti-gliadin-slgA antibodies of the samples are quantitatively determined by addition of a peroxidase-labeled anti-gliadin-slgA antibody mix. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. The enzymatic reaction is stopped by an acidic stop solution. The absorbance of the color compound is determined photometrically at 450 nm. The measured absorbance is directly proportional to the amount of bound anti-gliadin-slgA antibodies.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 9311MTP	PLATE	One holder with precoated strips	12 x 8wells
K 9311WB	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 9311K	CONJ	Conjugate (Peroxidase-labeled), ready to use	1 x 15 ml
K 9311KO1	CTRLNEG	Control negative, lyophilized	4 vials
K 9311KO2	CTRLPOS	Control positive, lyophilized	4 vials
K 9311CK	CTRLCUTOFF	Cut-off Control, lyophilized	4 vials
K 9311TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 9311AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

*Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0,2 µm) with an electrical conductivity of 0,055µS/cm at 25°C (≤18,2MΩ cm).

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be dissolved at 37°C using a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** could be stored in a closed flask at **2-8°C for one month**.

The diluted wash buffer is used as extraction buffer (see Sample preparation).

- The lyophilized **CTRLNEG, CTRLPOS and CTRLCUTOFF** (controls, negative, positive and cut-off) are stable at **2-8 °C** until the expiry date stated on the label. Reconstitution details are given in the data sheet. Diluted controls are not stable and can not be stored.
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SAMPLE PREPARATION

Extraction of the stool sample

1a. Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instruction for use

Please note that the dilution factor of the final stool suspension depends on the used amount of stool sample and the volume of the buffer.

SAS with 0.75 ml Buffer:

Applied amount of stool:	15 mg
Buffer Volume:	0.75 ml
Dilution Factor:	1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw Stool Sample has to be thawed. For remarkably inhomogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) **Fill the empty sample tube** with **0.75 ml** of ready-to-use extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into sample. The lower part of the dipstick exhibits notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off and leave 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for app. 10 minutes improves the result.
- e) Allow sample to stand for app. 10 minutes until sediment has settled down. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure, the sediment will not be dispersed again.

1b. Sample preparation kit from Roche Diagnostics, Mannheim, Germany (Cat. No. 10 745 804 322)

Alternatively, other stool sample preparation kits (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany) can be used. In the Roche sample preparation kit, 100 mg of stool sample are suspended in 5 ml of extraction buffer using a vibrator mixer (e.g. Vortex mixer). Centrifugation of the suspension is recommended.

Dilution Factor:**1:50**

100 µl of the resulting supernatant (**1a. or 1b.**) is used for the assay.

The stool suspension is not stable and can not be stored.

Raw stool can be stored at –20°C for 4 weeks. Avoid thawing and freezing cycles.

9. ASSAY PROCEDURE*Procedural notes*

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well.

Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.

Wash the PLATE (precoated microtiter plate) 5 x with 250 µl ELISA diluted wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.

Carry out the tests in duplicate.

1. Add **100 µl** of **CTRLNEG, CTRLPOS and CTRLCUTOFF** (controls, negative, positive and cut-off) and diluted **SAMPLE** (patient samples).
2. Incubate for **2 hours** at **2-8 °C**.
3. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** of CONJ (conjugate).
5. Incubate for **1 hour** at **2-8 °C**.
6. Decant the content of the plate and wash the wells **5 x with 250µl** ELISA wash buffer.
7. Add **100 µl** SUB (TMB substrate).
8. Incubate for **15 - 25 minutes** at room temperature.*
9. Add **50 µl** STOP (ELISA stop solution) and mix shortly.
10. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

10. RESULTS

Calculation

$$\text{Cut off} = 100 \text{ U/I} = OD_{\text{Cut off control}}^*$$

***The labeled parameters may need to be changed when running the test on automated systems. Please contact the producer for details.**

For example: cut off

Samples with a mean extinction higher than the OD of the cut off control are positive.

For example : positive patient

Patient extinction = 0.685

Extinction_{Cut off control} = 0.321 = 100 U/I

$$X = \frac{0.685 * 100}{0.321} = 213.4 \text{ U/I}$$

The calculation is for stool dilution of 1:50.

11. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

anti-gliadin-slgA-antibodies (stool): < 100 U/I

12. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) of the Immundiagnostik anti-gliadin-slgA-antibody-ELISA test was calculated from 20 replicate determinations on each one of two samples.

Intra-Assay CV n= 20

Sample	anti-gliadin-slgA-Ab mean value [U/l]	Intra-Assay CV [%]
1	527.7	2.9
2	348.9	3.9

The reproducibility of the results of two samples on different days was evaluated. The two samples were measured by different technicians on different days using the anti-gliadin-slgA-antibody-ELISA test.

Inter-Assay CV n= 12

Sample	anti-gliadin-slgA-Ab mean value [U/l]	Inter-Assay CV [%]
1	425.3	14.9
2	237.6	8.7

13. REFERENCES

1. Dieterich et al.: 1997; Nat. Med. 7 (3), 797
2. Rieken et al.: 1998; Dtsch. med. Wschr. 123, 1454
3. Green et al.: 1998; Clin. Persp. Gastroenterol. November, 133
4. Mothes, Th.: 1997; Münsch. Med. Wschr. 139, 111

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number