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anti-htTG slgA ELISA

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

*For the in vitro determination of anti-human tissue
transglutaminase slgA antibodies in stool*

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Additional information about our products is available on our homepage

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Tabel of content	Page
1. INTENDED USE	17
2. SUMMARY AND EXPLANATION OF THE TEST	17
4. MATERIAL SUPPLIED	19
5. MATERIAL REQUIRED BUT NOT SUPPLIED	19
6. PREPARATION AND STORAGE OF REAGENTS	20
7. PRECAUTIONS	20
8. SPECIMEN COLLECTION AND PREPARATION	21
9. ASSAY PROCEDURE	22
PROCEDURAL NOTES	22
TEST PROCEDURE	23
10. RESULTS	24
11. QUALITY CONTROL	24
EXPECTED VALUES	24
12. PERFORMANCE CHARACTERISTICS	25
PRECISION AND REPRODUCIBILITY	25
13. REFERENCES	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26

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Manual



anti-htTG sIgA ELISA

*For the in vitro determination of anti-human tissue
transglutaminase sIgA antibodies in stool*

Valid from 25.03.2011

K 9393

Immundiagnostik AG

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **anti-human tissue Transglutaminase slgA antibodies** in stool. For *in vitro* diagnostic use only.

2. INTRODUCTION

Dieterich et al. have shown, that tissue **transglutaminase** (tTG) is the predominant endomysial auto antigen characteristic for celiac disease.

In most celiac patients, usually several disorders are observed, e.g. malabsorption, infertility, osteoporosis and delayed growth of children. But it has also been widely reported that celiac disease is associated with a whole series of autoimmune diseases like Dermatitis herpetiformis Duhring, Diabetes mellitus, rheumatoid Arthritis, IgA-nephritis, neuro-psychiatric disorders, Hashimoto-Thyreoditis / M. Basedow and an increased risk of developing malignant T cell lymphoma.

Because the prevalence of associated autoimmune diseases in most cases is high, it is advisable to determine the auto antibodies against tissue **transglutaminase** (tTG) as a marker for celiac disease.


Indication

- Autoimmune disease
- Food intolerance

See also our proposal for stepwise diagnosis of gluten sensitivity on page 18.

3. PRINCIPLE OF THE TEST

The **Enzyme-Linked-Immuno-Sorbent-Assay** (ELISA) is used for semi-quantitative determination of anti-htTG slgA antibodies in faeces. In a first incubation step, the anti-htTG slgA antibodies in the sample are bound to their antigen (human, recombinant Transglutaminase), which is immobilized to the surface of the microtiter plates. To remove all unbound foreign substances, a washing step is carried out. In a second incubation step, a POD-labeled slgA antibody mix is added. After another washing step, to remove all unbound antibodies, the solid phase is incubated with the substrate, tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the amount of bound antibodies and can be determined photometrically at 450 nm.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 9393MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 9393WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 9393K	CONJ	Conjugate (Peroxidase-labeled), ready to use	1 x 15 ml
K 9393KO1	CTRLNEG	Control negative, lyophilized	4 vials
K 9393KO2	CTRLPOS	Control positive, lyophilized	4 vials
K 9393CK	CTRLCUTOFF	Cut-off Control, lyophilized	4 vials
K 9393TMB	SUB	TMB substrate (Tetramethylbenzidine)	1 x 15 ml
K 9393AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

*Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25°C (≤18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be dissolved at 37°C using a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** could be stored in a closed flask at **2-8°C for one month.**
The diluted wash buffer is used as extraction buffer (see Sample preparation, Point 8).
- The lyophilized **CTRLNEG, CTRLPOS and CTRLCUTOFF** (controls, negative, positive and cut-off) are stable at **2-8 °C** until the expiry date stated on the label. Reconstitution details are given in the data sheet. Diluted controls are not stable and cannot be stored.
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SAMPLE PREPARATION

Extraction of the stool sample

1a. Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instruction for use

Please note that the dilution factor of the final stool suspension depends on the used amount of stool sample and the volume of the buffer.

SAS with 0.75 ml Buffer:

Applied amount of stool:	15 mg
Buffer Volume:	0.75 ml
Dilution Factor:	1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw Stool Sample has to be thawed. For remarkably inhomogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) **Fill the empty sample tube** with **0.75 ml** of ready-to-use extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into sample. The lower part of the dipstick exhibits notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off and leave 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for app. 10 minutes improves the result.
- e) Allow sample to stand for app. 10 minutes until sediment has settled down. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure, the sediment will not be dispersed again.

1b. Sample preparation kit from Roche Diagnostics, Mannheim, Germany (Cat. No. 10 745 804 322)

Alternatively, other stool sample preparation kits (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany) can be used. In the Roche sample preparation kit, 100 mg of stool sample are suspended in 5 ml of extraction buffer using a vibrator mixer (e.g. Vortex mixer). Centrifugation of the suspension is recommended.

Dilution Factor: **1:50**

100 µl of the resulting supernatant (**1a. or 1b.**) is used for the assay.

The stool suspension is not stable and cannot be stored.

Raw stool can be stored at –20°C for 4 weeks. Avoid thawing and freezing cycles.

9. ASSAY PROCEDURE*Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any variations of the test procedure, without consulting the producer, may influence the test results. Immundiagnostik AG can therefore not be held responsible for any damage.
- The assay should always be performed according the enclosed manual.

Test procedure

Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well.

Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.

Wash the PLATE (precoated microtiter plate) 5 x with 250 µl ELISA diluted wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.

Carry out the tests in duplicate.

1. Add **100 µl** of **CTRLNEG, CTRLPOS and CTRLCUTOFF** (controls, negative, positive and cut-off) and diluted **SAMPLE** (patient samples).
2. Incubate for **2 hours** at **2-8 °C**.
3. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** of CONJ (conjugate).
5. Incubate for **1 hour** at **2-8 °C**.
6. Decant the content of the plate and wash the wells **5 x with 250µl** ELISA wash buffer.
7. Add **100 µl** SUB (TMB substrate solution).
8. Incubate for **15 - 25 minutes** at room temperature.*
9. Add **50 µl** STOP (stop solution) and mix shortly.
10. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

10. RESULTS

Calculation

$$\text{Cut off} = 100 \text{ U/l} = OD_{\text{Cut off control}}^*$$

***The labeled parameters may need to be changed when running the test on automated systems. Please contact the producer for details.**

For example: cut off

Samples with a mean extinction higher than the OD of the cut off control are positive.

For example : positive patient

Patient extinction = 0.685

Extinction_{Cut off control} = 0.234 = 100 U/l

$$X = \frac{0.685 * 100}{0.234} = 292.7 \text{ U/l}$$

The calculation is for stool dilution of 1:50.

11. QUALITY CONTROL

Control samples or stool pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of control samples, should be evaluated for plausibility using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable range, the results for the patient sample may not be valid.

Expected values

htTG-sIgA-antibodies in stool <100 U/l

12. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) of the Immundiagnostik anti-ht Transglutaminase-slgA-antibody-ELISA test was calculated from 20 replicate determinations on each one of two samples.

Intra-Assay CV n= 20

Sample	htTG-slgA Ab mean value [U/l]	Intra-Assay CV [%]
1	160.9	5.7
2	83.4	9.8

The reproducibility of the results of two samples on different days was tested. The two samples were measured by different technicians on different days using the anti-ht Transglutaminase-slgA-antibody-ELISA test.

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14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- All reagents in the kit package are for research use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Guidelines for medical laboratories should be observed.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number