anti-htTG IgA ELISA

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
For the in vitro determination of anti-human -tissue-Transglutaminase IgA antibodies in serum and plasma

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1. INTENDED USE	15
2. PRINCIPLE OF THE TEST	15
3. MATERIAL SUPPLIED	16
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. PREPARATION AND STORAGE OF REAGENTS	18
6. PRECAUTIONS	19
7. SPECIMEN COLLECTION AND PREPARATION	19
8. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
TEST PROCEDURE	20
9. RESULTS	21
10. LIMITATIONS	21
11. QUALITY CONTROL	21
EXPECTED VALUES	21
12. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
13. TECHNICAL HINTS	23
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Hersteller

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
anti-htTG IgA ELISA

For the in vitro determination of anti-human-tissue-Transglutaminase IgA antibodies in serum and plasma

Valid from 20.05.2011

K 9399

Immundiagnostik AG

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the determination of human **anti-tissue-Transglutaminase-IgA antibodies** in serum and plasma. For *in vitro* diagnostic use only.

2. PRINCIPLE OF THE TEST

This Enzyme Immuno Assay is a sandwich assay for the quantitative determination of anti-tissue transglutaminase antibodies (IgA anti-htTG). The wells of the microtiter plate are coated with the antigen. In a first incubation step, the anti-htTG IgA antibodies are bound to the coated antigen. To remove all unbound substances, a washing step is carried out.

In a further incubation step, a peroxidase-labeled antibody is added. After another washing step to remove all unbound substances, the solid phase is incubated with the peroxidase substrate, tetramethylbenzidine (TMB). Finally, an acidic stop solution is added to terminate the enzymatic reaction, whereby the colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the IgA anti-tissue-transglutaminase antibody concentration. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. IgA anti-tissue-transglutaminase antibodies present in the samples are determined directly from this curve.

Indication

- See our proposal for stepwise diagnosis of gluten sensitivity on page 17.

3. MATERIAL SUPPLIED


Cat. No	Content	Kit Components	Quantity
K 9399MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 9399WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 9399ST	STD	htTG-standards lyophilized	4 x 1 vial
K 9399K01	CTRL	Control, lyophilized	4 x 1 vial
K 9399K02	CTRL	Control, lyophilized	4 x 1 vial
K 9399K	CONJ	Conjugate (rabbit-anti-IgA antibody, peroxidase labeled)	1 x 200 µl
K 9399TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 9399AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker with 37 °C incubator
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

*Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0,2 µm) with an electrical conductivity of 0,055µS/cm at 25°C (≤18,2MΩ cm).

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5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be re-dissolved at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **CONJ** (conjugate) must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The undiluted conjugate is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- The lyophilized **STD** (standard) and **CTRL** (controls) must be reconstituted with **500 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes at room temperature, and mix thoroughly by gentle inversion to insure complete reconstitution. The undiluted standard and controls are stable at **2-8 °C** until expiry date stated on the label. **Diluted standard and controls are not stable and cannot be stored.**

The preparation of the **standard curve** (volumes and concentrations) is described in the product specification.

- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

6. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, all kit components should be handled as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

7. SPECIMEN COLLECTION AND PREPARATION

Serum/Plasma

Dilute serum and plasma **1:250** in ELISA wash buffer in two dilution steps.

Example:

Dilution step 1: **10 µl serum + 90 µl wash buffer**

Dilution step 2: **40 µl dilution 1 + 960 µl wash buffer**

8. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Perform assay always according the enclosed manual.

Test procedure

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** of **STD** (standards), **CTRL** (controls) and **diluted samples** per well.
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the plate content and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** of diluted **CONJ** (conjugate, peroxidase labeled anti-IgA antibody) into each well.
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the plate content and wash the wells **5 x with 250 µl** ELISA wash buffer.
7. Add **100 µl SUB** (TMB substrate solution).
8. Incubate for **10 - 20 minutes** at room temperature.
9. Add **50 µl STOP** (stop solution) and mix shortly.
10. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

9. RESULTS

Calculation

Since the sample dilution of 1:250 is already considered in the calibration curve, the dilution factor is 1.

If a sample dilution different from 1:250 has been used, a corresponding dilution factor must be considered, e.g., at a dilution of 1000 the results must be multiplied by 4.

10. LIMITATIONS

Samples with IgA anti-htTG levels higher than the highest standard should be further diluted with wash buffer and re-assayed.

11. QUALITY CONTROL

Controls should be analyzed with each run. Results, generated from the analysis of controls, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality controls are outside the acceptable limits.

Expected values

htTG IgA Serum cut off > 22 AU/ml positive

htTG IgA Serum cut off < 16 AU/ml negative

Results between 16 - 22 AU/ml are in the grey range.

We recommend each laboratory to establish its own normal range.

In order to avoid unnecessary biopsies, tTG-IgA titers within the grey zone should be controlled after 3 - 6 months. Biopsies should be performed only in case of increased tTG-IgA titers. Until now, there is no information on correlation between age and titer level in children with negative tTG-IgA values. Age-normalized threshold values could be useful.*

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Precision and reproducibility

Intra-Assay-Variation

The intra-assay variation of the assay was calculated from 64 determinations of one sample by one person.

Intra-Assay CV n = 64

Sample	IgA anti-htTG mean value [AU/ml]	Intra-Assay CV [%]
1	34,12	5,2

Inter-Assay-Variation

The inter-assay variation of the assay was calculated from data on 2 samples obtained in 15 or 17 determinations by different persons on two different days.

Inter-Assay CV n = 15/17

Sample	IgA anti-htTG mean value [AU/ml]	Inter-Assay CV [%]
1 (n=15)	15,74	13,7
2 (n=17)	68,04	15,4

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in-vitro* diagnostics only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number