
MutaREX[®] Norovirus *real time* RT-PCR Kit

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60°C für **30 sec**

72°C für **20 sec**

ramping time: **20°C/ sec** – acquisition mode here:

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

MutaREX[®] Norovirus

real time RT-PCR Kit

Qualitative assay for specific detection of **Norovirus** (GG I, II) using *real time* PCR capillary systems (e. g. LightCycler[®] 1.5/ 2.0, **Roche***) and parallel control of RNA extraction efficiency from stool samples.

REF KG290132 

REF KG290196 

* MutaREX[®] Norovirus is **licensed** from Roche Molecular Systems, Inc.

For in vitro Diagnostic use only

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany
www.immundiagnostik.com
info@immundiagnostik.com

Tel.: +49 (0)6251/ 701900
 Fax: +49 (0)6251/ 849430

1. INTENDED USE

The **MutaREX[®] Norovirus** *real time* RT-PCR kit is a qualitative assay for the specific detection of Norovirus (GG I and II) in stool samples using *real time* RT-PCR capillary systems (e.g. LightCycler[®] 1.5/ 2.0, Roche). In parallel, the RNA extraction efficiency from stool samples is controlled by a separate internal control. The kit is **licensed** from **Roche Molecular Systems, Inc.**

2. INTRODUCTION

Non-bacterial **Gastroenteritis** is mainly caused by species of human-pathogen noroviruses. These high infectious viruses have a non-enveloped RNA genome and belong to the family of **Caliciviridae**. They are resistant against environmental influences, e.g. against temperatures between -20°C and +60°C, acids (pH3) and chlorite (10mg/ L). NV infections may occur the whole year but in Europe a seasonal increase from October - March is typical. NV trans-mission happens **direct** from person to person by contact-/ smear-infection or also **indirect** by contaminated food or water. Already a low dose of NV particles is sufficient to release epidemic infections. Therefore especially closed environments (hospitals, elderly homes, kindergartens, sport-courts, cruise-ships etc.) are affected. Positive tested persons have to be registered (RKI) and isolated for several days.

Due to their high genetic variability, human-pathogen NV is subdivided in different geno-groups/ genotypes. But because **MutaREX[®] Norovirus** reagents are routine up-dated using different sequence data bases, they recognize even strains from high genetic diversity, e.g.:

GGI: Norwalk, Southampton, Queens Arms, Desert Shield, Winchester, Sindlesham, Chiba

GGII: Lordsdale, Melksham, Hawaii, Mexico, Leeds, Hillingdon, Snow Mountain, Toronto, Bristol

3. PRINCIPLE OF THE TEST

MutaREX[®] Norovirus *real time* RT-PCR kit contains specific primers, fluorescence-marked probes and additional material for the detection of Norovirus (genogroup I and II). Human stool samples may be used as starting material for RNA-extraction and following pathogen analysis with **MutaLREX[®] Norovirus**.

Reverse transcription (RT) of potentially contained viral RNA to cDNA and subsequent amplification of norovirus-specific fragments by means of PCR is done in only **one step**. **Target sequence** for the detection is the region in the **ORF1 / ORF2 - junction**.

Detection of NV amplicates is achieved in *real time* by hybridization and hydrolysis of NV-specific fluorescence probes. Their emitted signal is measured (during PCR process) by the optical unit of the *real time* PCR system in use (LightCycler 1.5 respectively 2.0- software). **Norovirus-specific amplification** is measured in the **F1**-channel (**530** nm).

Furthermore, **MutaREX[®] Norovirus** *real time* RT-PCR - kit contains also a separate RNA extraction control. This is added initially prior the nucleic acid extraction procedure and detected in the follow-up by an independent (heterologous) amplification step. This internal control (IC) allows on the one hand the exposure of errors during RNA extraction procedure and identifies on the other hand a potential inhibition of reverse transcription or PCR. Doing this, the risk for false-negative results is minimized extensively. **RNA extraction control** (= internal control for PCR) is measured at **710** nm in channel **F3** (LightCycler 1.5) respectively **F6** (LightCycler 2.0).

4. KIT CONTENT

Each kit contains enough reagents to perform 32 or 96 tests and also a package insert.

	Reagent	Confection 32		Confection 96		Colour
A1	Enzyme Mix	1x	30 µl	1x	90 µl	blue
A2	Primer-/ Probe Mix	1 x	500 µl	2 x	750 µl	yellow
A3	Positive Control	1 x	20 µl	2 x	30 µl	red
A4	Negative Control	1 x	100 µl	1 x	200 µl	green
A5	Extraction Control (IC)	2 x	90 µl	3 x	200 µl	transparent

5. TEST PERFORMANCE

Required materials - provided:

- Reagents for *real time* RT-PCR
- Package insert

Required materials - not provided:

- *real time* PCR capillary system (e. g. LightCycler® instrument, Roche)
- *real time* PCR reaction tubes (e. g. LightCycler® capillaries, Roche)
- Table centrifuge (e. g. LightCycler® capillary centrifuge, Roche)
- Cryocontainer for PCR reaction tubes (e. g. LightCycler® cooling block, Roche)
- Color Compensation Kit for the used *real time* PCR system
- RNA extraction kit for stool samples (e. g. **MutaCLEAN Stool (Viral RNA)**, KG1035, Immundiagnostik)
- Pipets (0.5 µl – 1000 µl) with sterile filter tips
- sterile reaction tubes (1,5 ml and 2,0 ml)

6. STORAGE AND HANDLING

- All reagents (A1 to A5) should be **stored till immediate use** at **<-20°C**.
- **Avoid several freeze / thawing** cycles for the reagents A2, A3 and A5 (if case of sporadic use prepare **suited aliquots**).
- **cool all reagents** during the working steps.
- Primer-/ Probe-Mix (A2) should be stored dark (**light protection of** fluorescence probes).
- All reagents can be used until the expiration date printed on the labels.

7. WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- **MutaREX® Norovirus** is **Roche-licensed**, no further fees for „reported-results“ necessary.
- This assay needs to be carried out by in Molecular Biology methods skilled personnel.
- Clinical samples should be regarded as potentially infectious materials.
- This assay needs to be run according to GLP (Good Laboratory Practice).
- Do not use the kit after its expiration date.

AMPLIFICATION: The PCR technology is utmost sensitive. Thus, amplification of a single molecule generates millions of identical copies. These copies may evade through aerosols and sit on surfaces. In order to avoid contamination of samples with RNA which previously was amplified, it is important to physically strictly divide sample and reagent preparation units from sample amplification units. Set up two separate working areas:

- 1) Aerea for sample preparation.
- 2) Aerea for preparation of Master Mix.
- 3) Aerea for pipetting Master Mix and samples to the PCR tubes.

Pipets, vials and other working materials should not circulate among working units!

- Use always sterile pipette tips with filters and avoid aerosols.
- Wear separate coats and gloves in each area.
- Routinely decontaminate your pipettes and laboratory benches with decontaminant (do not use ethanol solutions).

8. PROCEDURE

The complete procedure is separated in three steps:

- A) RNA extraction from stool samples (e.g. with **MutaCLEAN® Stool (Viral RNA)**, KG1035 from **Immundiagnostik**).
- B) reverse transcription of RNA and subsequent amplification/ combined detection of generated cDNA templates using hybridisation probes.
- C) Result interpretation with the software of *real time* PCR instrument (LightCycler 1.5/ 2.0).

A) RNA EXTRACTION

- 1) extraction of viral RNA is done with a **commercial** available RNA isolation kit (e. g. **MutaCLEAN® Stool (Viral RNA)**, KG1035 from **Immundiagnostik**) from clinical stool samples according to the manufacturer`s instructions.

For doing this, resuspend a “pea-sized” stool sample in 1,5 ml **pure water** (do **not** use buffers!) by intensive vortexing and let sediment solid particles (if necessary support by centrifugation for 2 min at 5000 rpm in table centrifuge). Use the supernatant for RNA extraction without any further dilution. Very liquid stool samples should be used directly (without further resuspension) for the RNA extraction, but also here vortex intensively and proceed as before described.

- 2) Important: **Additional** to the samples should be extracted also a water control which is for identification of eventually occurrence of inhibiting matrix effects. Treat this water control analogous to the samples.
- 3) Extraction control (A5) in the **MutaPLEX® Norovirus** – Kit is added to each sample BEFORE begin of the extraction procedure: Prepare a “working solution“ with the first buffer of the used RNA isolation kit (eventually plus supplements as „Poly A“) and the extraction control (A5). Multiplicate the buffer volume necessary for one extraction (this is 200 µl for **MutaCLEAN® Stool (Viral RNA)** from **Immundiagnostik**) with number (n) of extractions performed in total (do not forget the water control and add also one virtual extra volume for reasons of unprecise pipetting) and add for each extraction 5µl of A5:
= [extraction buffer volume + 5µl A5] x [n+1]. Perform with that “working solution” the extraction procedure according to the manufacturer`s manual.

OPTIONAL: If no control of extraction efficiency is needed, the reagent A5 may be used also directly in the PCR Master Mix as internal control (IC). In that case use 0.5 µl per reaction!

- 4) If *real time* RT-PCR is not performed immediately, store the extracted RNA according to the manufacturer's recommendation (normally at < -20°C).

B) *Real time* Norovirus (MutaREX®) RT-PCR-PROTOCOL

Please **read** carefully the manufacturer's instructions **before** starting the procedure! Each assay should include a negative and positive control. Use filter tips for all pipetting steps. The necessary total Master Mix volume is calculated as follows:

- 1) Multiply Enzyme Mix volume per reaction with the number (n) of necessary reactions (samples and controls) + add one extra virtual volume for reasons of unprecise pipetting; proceed in the same manner with all additional reagents!

Cool all reagents during the working steps!

- 2) For **PCR Master Mix** preparation combine the following reagents in a sterile tube and mix gently by pipetting (about **15 - 20x !** - do **NOT vortex!**):

Reaction Volume	Master Mix Volume
0.8 µl Enzyme Mix (A1)	0.8 µl x (n+1)
14.2 µl Primer-/ Probe Mix (A2)	14.2 µl x (n+1)

- 3) Pipet **15 µl** of this Master Mix using micropipets with sterile filter tips in each of the *real time* PCR reaction tubes (e. g. LightCycler® capillaries); do not forget the controls.

Add **5 µl** of sample RNA (or positive/ negative/ water control) to each corresponding PCR reaction tube (LightCycler capillary). It is recommended to close the capillaries immediately after filling and to prepare the negative control indeed first but to close at last (as contamination control).

Spin briefly to collect all reagents at the bottom of capillaries (use LC capillary centrifuge).

- 4) Transfer the capillaries into the *real time* PCR system (LightCycler rotor) and run the RT-PCR using following temperature protocol:

50°C for **10 min** (reverse transcription)

95°C for **10 min**

40 cycles: (amplification)

95°C for **15 sec**

60°C for **30 sec**

72°C for **20 sec**

40°C for **10 sec**

ramping time: **20°C/sec** – acquisition mode here:

NONE

SINGLE

NONE


C) PCR ANALYSIS AND INTERPRETATION OF RESULTS

- 1) Perform *real time* RT-PCR (in LightCycler® capillary system, Roche) with **MutaREX® Norovirus**.
- 2) Switch **on** the **color compensation** filter required because of the simultaneous use of several fluorescence-marked probes (by activating the field "*Choose CCC File*").

- 3) The result of **specific Norovirus amplification** is shown in the **F1-channel** (at **530 nm**) and the result of the extraction efficiency control (= internal control **IC**) is shown at **705 nm** (channel F3 respectively F6).

Following results can arise:

- A signal is detected at 530 nm (F1).
The result is positive: **The sample contains Norovirus RNA.**
In this case, the detection of a signal at 705 nm is inessential, as high concentrations of *Norovirus* cDNA can lead to a reduced or absent fluorescence signal of the internal control in channel F3 (competition).
- No signal is detected at 530 nm (F1), but only at 705 nm (F3, signal of the internal control).
The result is negative: **The sample does not contain any Norovirus RNA.**
The detected signal of the internal control excludes the possibility of an inhibition.
- Neither at 530 nm (F1) nor at 705 (F3) nm is a signal detected.
A diagnostic statement can not be made.
An **inhibition** of the *real time* RT-PCR occurred.

Note:

- PC is concentrated with about 100.000 copies/ μ l (means 500.000 in PCR) and generates amplification signals with ct-value 24 – 25 (+/- 1). Each 1:10 – dilution results in a (logarithmic) ct-value increase of 3.3 units.
- Sensitivity of the PCR is about 1 – 5 copies/ μ l.