



**LDN Labor Diagnostika Nord GmbH & Co. KG**

Am Eichenhain 1, 48531 Nordhorn

Telefon: +49-5921-8197 0

Telefax: +49-5921-8197 222

e-mail: [info@ldn.de](mailto:info@ldn.de)

Internet: <http://www.ldn.de>

**LDN<sup>®</sup>**

## Instructions for use **a-Amylase Saliva**

**REF**

**SA P-6900**



**RUO**

For Research use only-  
Not for use in diagnostic  
procedures

## **INTENDED USE**

The  $\alpha$ -Amylase Saliva assay is a kinetic colorimetric method for quantitative determination of  $\alpha$ -amylase in saliva.

## **Clinical Significance**

Amylase is the name given to enzymes that break down starch. They are classified as saccharidases, enzymes that cleave polysaccharides.

Although the amylases are designated by different greek letters, they all act on  $\alpha$ -1,4 glycosidic bonds.

The  $\alpha$ -amylases are calcium metalloenzymes, completely unable to function in the absence of calcium. The  $\alpha$ -amylase breaks down long-chain carbohydrates, ultimately yielding maltotriose and maltose from amylose, or maltose, glucose and "limit dextrin" from amylopectin. In animals, it is a major digestive enzyme.

Although found in many tissues, the  $\alpha$ -amylase is most prominent in pancreatic juice and saliva which each have their own isoform.

Salivary  $\alpha$ -amylase breaks starch down into maltose and dextrin. This form is also called ptyalin. Ptyalin will break large, insoluble starch molecules into soluble starches (amylodextrin, erythrodextrin, achrodextrin) producing successively smaller starches and ultimately maltose. Ptyalin acts on linear  $\alpha$ (1,4) glucosidic linkages, but compound hydrolysis requires an enzyme which acts on branched products. Salivary amylase is inactivated in the stomach by gastric acid.

Pancreatic  $\alpha$ -amylase randomly cleaves the  $\alpha$ (1-4)glycosidic linkages of amylose to yield dextrin, maltose or glucose molecules.

## **PRINCIPLE**

The human  $\alpha$ -Amylase hydrolyses the 2-chloro-4 nitrophenyl- $\alpha$ -maltotrioside (CNP-G3) in glucose polymers and short-chain p-nitrophenyl-oligosaccharide with formation of 2-chloro-4-nitrophenol (CNP).

The increase of the extinction is evaluated spectrophotometrically at 405 nm and is proportional to  $\alpha$ -amylase activity in the sample.

## **Reagent, material and instrumentation**

### **Reagent and material supplied in the kit**

**RELEASE-BUFF**      **SA P-6926 Reagent A (Release Buffer)**

(1 bottle) 32 mL; CNP-G3 2 mmol/L, Goods buffer 100 mmol/L, stabilisers and preservatives

**ASSAY-BUFF**      **SA P-6913 Conc. Assay Buffer 5x**

(1 bottle) 50 mL; Hepes buffer 200 mM pH 7.4; BSA 0,5 g/L

**96**      **BA D-0032 Microplate**

### **Reagents necessary not supplied**

Distilled water.

### **Auxiliary materials and instrumentation**

Automatic dispenser.

Microplates reader (405 nm).

**Note:** Store all reagents between 2°C – 8°C in the dark.

## **PRECAUTION**

- Once open, store all reagents between 2°C – 8°C and do not use them beyond the expiration date.
- For in vitro diagnostic use only.
- Avoid the contact with reagents which could be toxic if are ingested. Not pipette with the mouth.

## **PROCEDURE**

### **Preparation of Assay Buffer**

Dilute the whole bottle of *Assay Buffer Conc 5x* in 200 mL of distilled or deionized water.  
Keep between 2°C – 8°C until expiration date (see vial label).

### **Preparation of the Sample**

For sample collection it is advised to use glass centrifuge tubes and plastic straws.  
Do not use plastic tube or commercially available devices for the saliva collection to avoid false results.  
Let the saliva flow down through the straw into the centrifuge glass tube, freeze and thaw it to help to mucins precipitation. Centrifuge at 3000 rpm for 15 minutes.

Dilute 10 µL of liquid supernatant to 1 mL of diluted Assay Buffer (reagent 2).  
Mix gently by leaving it for at least 5 minutes on a rotating shaker.

If the assay is not carried out in the same day of collection, store samples at -20°C.

## **PROCEDURE**

All the reagents should be brought to room temperature 22°C - 28°C.

If you should manually dispense a high number of samples, it is advised to use maximum four strips for each tests.

Format the microplate wells for each patient specimen to be assayed in duplicate. Replace any unused microwell strips back into the aluminium bag, seal off and store at 2°C – 8°C

	Sample	Blank
Diluted Sample	10 µL	---
Distilled water	---	10 µL
Reagent A	300 µL	300 µL

Incubate at 37°C for 3 min immediately after reagents dispensation.

At the end of incubation time put the microplate on microplate reader at ambient temperature (22-28°C) and read the **absorbance variation (ΔA)** two times at 405 nm, the first after 1 minute and the second after 5 minute from the end of incubation time, subtracting each time the absorbance of blank.

## **QUALITY CONTROL**

Each laboratory should assay controls at normal, high and low levels range of α-amylase for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## **LIMITATION OF PROCEDURE**

### **Assay Performance**

It is important that the time of reaction in each well is held constant for reproducible results. Plate readers measure the OD vertically. Do not touch the bottom of the wells.

## **RESULTS**

Calculate  $\Delta A/\text{Min.} = [(A_{5 \text{ min}} - A_{1 \text{ min}}) / 4]$

Calculate:  $\alpha$ -Amylase (U/mL):

$$\frac{\Delta A/\text{min}}{\text{MMA}} \times \frac{\text{TV}}{\text{SV}} \times \frac{\text{DF}}{\text{LP}}$$

Where :

$\Delta A/\text{Min}$  = Absorbance difference per minute

TV = Total assay volume (0.310 mL)

DF = Dilution factor

MMA = millimolar absorption (extinction) coefficient of 2-chloro-p-nitrophenol (12.9)

SV = Sample volume (0.010 mL)

LP = Light path in centimeters = 0.95 (specific to microplate received with the kit)

$$\frac{\Delta A/\text{min}}{12,9} \times \frac{0,310}{0,010} \times \frac{100}{0.95} = \Delta A/\text{min} \times 253 = \text{U/ml} (*)$$

(\*) =  $\alpha$ -amylase activity in sample

### **Note:**

The microplate provided with the kit is not chemically treated and its use is intended like a batch of reaction. Otherwise the reaction may be conducted in a different batch (for example in a spectrophotometer cuvette) and it is important to know the right optical path length (LP) that is the distance the light travels through the sample; for standard cuvette it is generally 1 cm.

## **Reference Value**

Adults (Normal) 80.0 U/mL

Absolute Range: 4 – 400 U/mL

## **Performance and Characteristics**

### **Precision**

#### Intra-Assay

Within run variation was determined by replicate determination (16x) of three different levels of control sera in one assay. The within assay variability is < 1.5 %.

#### Inter-Assay

Between run variation was determined by replicate measurements of three different saliva samples in 2 different lots. The between assay variability is <1.5 %.

### **Sensitivity**

The lowest detectable concentration of  $\alpha$ -amylase that can be distinguished from the zero standard is 2.5 U/mL at the 95 % confidence limit.

### **Correlation**

The  $\alpha$ -Amylase kit (ENZ-5127) was compared to a similar commercially available kit.

The linear regression curve was calculated

$$y = 1.19 x + 2.55$$

$$r = 0.998 (r^2 = 0.996)$$

## **WASTE MANAGEMENT**

Reagents must be disposed off in accordance with local regulations.

## **REFERENCES / Literature**

1. Tietz NW. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia: WB Saunders Company: 46-49 (1995)
2. (IFCC) "Approved Recommendation on IFCC Methods for the measurement of catalytic concentration of Enzymes" \_Part 9
3. Clin Chem Lab Med 36(3): 185 –203 (1998)

## **TROUBLESHOOTING**

### **ERRORS / POSSIBLE CAUSES / SUGGESTIONS**

#### **No colorimetric reaction**

- no Reagent A pipetted
- contamination of Reagent A
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

#### **Too low reaction (too low ODs)**

- incorrect Reagent A (e.g. not from original kit)
- incubation time too short, incubation temperature too low

#### **Too high reaction (too high ODs)**

- incorrect Reagent A (e.g. not from original kit)
- incubation time too long, incubation temperature too high

#### **Unexplainable outliers**

- contamination of pipettes, tips or containers

#### too high within-run CV%













- reagents and/or strips not pre-warmed to room temperature prior to use
- incubation conditions not constant (time, temperature)

#### too high between-run CV %

- too long dispense time
- person-related variation

---

### **Symbols:**

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE labelled
	Caution		Catalogue number		For research use only!