



Sanquin

# PeliKine compact™ human IL-13

REF

M1913

## Introduction

Interleukin 13 (IL-13) is an immunoregulatory protein produced by activated T-cells. The protein encoded by the IL-13 cDNA is the human homologue of a mouse TH2 cytokine product called P600 [1-3]. IL-13 shares many of its biological activities with the TH2 cytokine Interleukin 4; both cytokines are able to enhance expression of CD23 on monocytes and B-cells and also induce IgE production [3-8]. Production of many LPS-induced monokines, such as IL-1 $\alpha$ , IL-1 $\beta$ , IL-6, IL-8, IL-10, TNF $\alpha$ , IFN $\alpha$ , MIP-1 $\alpha$ , GM-CSF and G-CSF are inhibited by IL-13 [3,6], whereas IL-1ra is upregulated. These properties are shared with IL-4 and IL-10. Therefore IL-4 and IL-13 can be considered as anti-inflammatory molecules. In contrast to IL-4, IL-13 has no growth-promoting effect on T-cells and cannot compete for IL-4 binding to a human T-cell line. Therefore it was thought that the specific receptor for IL-13 is lacking on T-cells [4,9]. However, recently an inhibitory effect of IL-13 on IL-8- and RANTES-induced chemotaxis of T-cells has been described, indicating that T-cells do respond to IL-13 [10], possibly by inhibition of production of the TH1 inducer IL-12 [4]. Bioassays for the quantification of IL-13, including the proliferative assay of an IL-13- dependent subclone of the B9 cell line [12] can be used. However, IL-13 shares many of the biological effects of IL-4 and for this reason the two commonly interfere in bioassays. Furthermore, these bioassay, although sensitive, are time consuming and susceptible to interference by other substances. The PeliKine compact™ human IL-13 ELISA kit [11] has been developed for faster, more reproducible and specific quantification of human IL-13 in serum, plasma and other body fluids, as well as in cell-culture supernatant.

## Assay procedure

See Assay procedure for PeliKine™ compact ELISA kit.

## Kitcomponent list

Quantity	Kit component	Volume	Cap colour
1 vial	coating antibody	100-fold concentrated	375 $\mu$ l
1 vial	blocking reagent	50-fold concentrated	2 ml
1 vials	IL-13 standard (lyophilized)	see label	500 $\mu$ l
1 vial	biotinylated antibody	100-fold concentrated	375 $\mu$ l
1 vial	streptavidin-poly-HRP conjugate	10,000-fold concentrated	20 $\mu$ l
1 bottle	HPE-dilution buffer	5-fold concentrated	55 ml
3 pcs	microtiter plate + lid	-	-
10 pcs	plate seals	-	-

## Sensitivity

MEAN calculated zero signal + 3 SD : 0.5 – 1.5 pg/ml (shake – static incubation)  
2x (MEAN calculated zero signal) : 1 – 3 pg/ml (shake – static incubation)

## Expected values

IL-13 values in fresh serum and plasma samples of healthy individuals are below 10 pg/ml.

## Specificity

No crossreactivity was observed with the following recombinant human proteins: IL-1 $\alpha$ , IL-1 $\beta$ , IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-9, IL-10, IL-11, Macrophage Colony Stimulating Factor (M-CSF), Granulocyte Colony Stimulating Factor (G-CSF), Granulocyte/Macrophage Colony Stimulating Factor (GM-CSF), Leukaemia Inhibitory Factor (LIF), RANTES, Stem Cell Factor/ Mast Cell Factor (SCF/MCF), Transforming Growth Factor  $\beta$ -1 (TGF $\beta$ -1), Tumour Necrosis Factor  $\alpha$ , Tumour Necrosis Factor  $\beta$  (TNF $\beta$ /Lymphotoxin), and Interferon  $\gamma$  (IFN $\gamma$ ).

## References

1. Mosmann, T.R. and Coffman, R.L. (1989) Annu.Rev.Immunol. **7**: 145
2. Minty, A. *et al* (1993) Nature **362**: 248
3. McKenzie, A.N.J. *et al* (1993) Proc.Natl.Acad.Sci. USA **90**: 3735
4. Zurawski, G. and De Vries, J.E. (1994) Immunol.Today **15**: 19
5. Zurawski, G. and De Vries, J.E. (1994) Stem Cells **12**: 169
6. De Waal Malefyt, R. *et al* (1993) J.Immunol. **151**: 6370
7. Defrance, T. *et al* (1994) J.Exp.Med. **179**: 135
8. Punnonen, J. *et al* (1993) Proc.Natl.Acad.Sci. USA **90**: 3730
9. Zurawski, G. *et al* (1993) EMBO J. **12**: 2663
10. Jinquan, T.B. *et al* (1995) J.Immunol. **154**: 3742
11. Van der Pouw Kraan, C.T.M. *et al* (1995) submitted
12. Labit-Le Bouteiller, C. *et al* (1995) J.Immunol.Meth. **181**: 29

Reagents  
Plesmanlaan 125  
1066 CX Amsterdam  
The Netherlands  
Phone: +31 20 512 3599  
Fax: +31 20 512 3570  
E-mail: reagents@sanquin.nl  
Website: [www.sanquinreagents.com](http://www.sanquinreagents.com)

p.t.o.  
M1913/01 1606090843

Specification sheet

**Standard** A natural human IL-13 standard has been calibrated against the WHO International Standard (IL-13 94/622; National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, U.K. 1 WHO Unit = 1 ng IL-13).

The kit contains one lyophilized vial with natural IL-13.


Reconstitute the lyophilized standard by adding 500 µl of distilled water to the vial. Incubate for 10 minutes at room temperature and mix gently. After reconstitution the standard must be kept cold (2-8°C) and stored frozen after use (<-18°C, preferably <-70°C).

**Standard curve** Label 7 tubes, one tube for each dilutions: 125, 50, 20, 8, 3.2, 1.3 and 0.5 pg/ml. Pipette 420 µl of working-strength dilution buffer into the tube labelled 125 pg/ml and 300 µl of workingstrength dilution buffer into the other tubes. Transfer 60 µl of the IL-13 standard (1000 pg/ml) into the first tube labelled 125 pg/ml, mix well and transfer 200 µl of this dilution into the second tube labelled 50 pg/ml. Repeat the serial dilutions five more times by adding 200 µl of the previous tube of diluted standard to the 300 µl of dilution buffer. The standard curve will contain 125, 50, 20, 8, 3.2, 1.3, 0.5 and 0 pg/ml (dilution buffer).

It is recommended to prepare two separate series for each assay.

**Samples** It is recommended to dilute the test samples at least 1:2 in working-strength dilution buffer. If high levels of IL-13 (outside the standard curve) are expected in the test samples, additional dilutions of sample i.e. 1:20 and 1:100 should also be prepared.

**Typical standard curve**



	STATIC INCUBATION	SHAKEN INCUBATION
	Calculated mean absorbance at 450 nm	
substrate blank	0	0
0 pg/ml	0.036	0.065
0.5 pg/ml	0.053	0.104
1.3 pg/ml	0.081	0.144
3.2 pg/ml	0.158	0.289
8 pg/ml	0.314	0.573
20 pg/ml	0.808	1.394
50 pg/ml	1.627	2.772
125 pg/ml	2.030	≥ 3.000

**DO NOT USE THESE DATA TO CONSTRUCT A STANDARD CURVE FOR SAMPLE VALUE CALCULATIONS**