

Sanquin

M1916

Introduction

REF

Interleukin 6 (IL-6) is a mediator of the inflammatory response and is involved in the induction of acute phase proteins [1,2,3,4] and the development of fever [5]. A marked correlation between IL-6 levels and inflammatory processes has been demonstrated in synovial fluid and serum of rheumatoid arthritis patients [6,7,8] and in serum of patients with burns [9,10]. It was demonstrated that in recipients of kidney transplants the IL-6 levels in serum and urine hallmark the onset of rejection episodes [11,12]. Elevated IL-6 levels were also observed in sera of patients with septic shock, multiple myeloma and alcoholic hepatitis, and a significant difference between IL-6 levels of survivors and non-survivors was found [13,14,15].

Bioassays for the quantification of IL-6, based on the proliferation of B-cell hybridomas have been used for several years [16,17,18]. These assays, although sensitive, are time consuming and susceptible to interference by other substances.

This PeliKine compact[™] IL-6 ELISA kit [19] has been developed for faster, more reproducible and specific quantification of human IL-6 (huIL-6) in plasma and other body fluids, as well as in cell-culture supernatant.

Assay procedure

See Assay procedure for PeliKine[™] compact ELISA kit.

Kitcomponent list

Quantity	Kit component		Volume	Cap colour
1 vial	coating antibody	100-fold concentrated	375 <i>µ</i> l	red
1 vial	blocking reagent	50-fold concentrated	2 ml	transparent
1 vials	IL-6 standard	see label	750µl	black
1 vial	biotinylated antibody	100-fold concentrated	375 <i>µ</i> l	yellow
1 vial	streptavidin-poly-HRP conjugate	10,000-fold concentrated	20 <i>µ</i> I	brown
1 bottle	HPE-dilution buffer	5-fold concentrated	55 ml	
3 pcs	microtiter plate + lid	-	-	
10 pcs	plate seals	-	-	

IL-6 values in fresh serum and plasma samples of healthy individuals are below 20 pg/ml.

No crossreactivity was observed with the following recombinant human proteins: IL-1α, IL-1β, IL-2, IL-3, IL-4, IL-5, IL-7, IL-8, IL-9, IL-10, IL-11, IL-13, Macrophage Colony Stimulating Factor (M-CSF), Granulocyte Colony Stimulating Factor (G-CSF), Granulocyte/Macrophage Colony Stimulating Factor (GM-CSF), Leukemia Inhibitory Factor (LIF), RANTES, Stem Cell Factor/ Mast Cell Factor (SCF/MCF), Transforming Growth Factor β-1 (TGFβ-1), Tumour Necrosis Factor α (TNF-α), Tumour

0.2 - 0.4 pg/ml (shake - static incubation)

0.5 - 1.0 pg/ml (shake - static incubation)

Expected values

Specificity

Sensitivity

References

1. Gauldie, J. et al (1987) Proc.Natl.Acad.Sci.(USA) 84: 7251

2. Le, J. and Vilcek, J. (1989) Lab.Invest. <u>61</u>: 588


Factor ß (TNFB/Lymphotoxin), and Interferon γ (IFN γ).


- 3. Heinrich, P.C. et al (1990) Biochem. J. 265: 621
- 4. Kishimoto, T. (1989) Blood 74: 1

MEAN calculated zero signal + 3 SD : 2x (MEAN calculated zero signal) :

- 5. Helle, M. et al (1988) Eur. J. Immunol. 18: 957
- 6. Houssiau, F.A. et al (1988) Arthritis Rheum. 31: 784
- 7. Swaak, A.J.G. et al (1988) Scand. J. Rheumatol. 17: 469
- 8. Waage, A. et al (1989) Clin.Immunol.Pathol. 50: 394
- 9. Nijsten, M.W.N. et al (1988) Lancet 11: 921
- 10. Guo, Y. et al (1990) Clin.Immunol.Pathol. 54: 361
- 11. Van Oers, M.H.J. et al (1988) Clin.Exp.Immunol. 71: 314
- 12. Yoshimura, N. et al (1991) Transplantation 51: 172
- 13. Hack,C.E. et al (1989) Blood 74: 1704
- 14. Ludwig, H. et al (1991) Blood 77: 2794
- 15. Sheron, N. et al (1991) Clin. Exp. Immunol. 84: 449
- 16. Aarden, L.A. et al (1985) Lymphokines 10: 175
- 17. Van Snick, J. et al (1987) J.Exp.Med 165: 641
- 18. Helle, M. et al (1988) Eur. J. Immunol. 18: 1535
- 19. Helle, M. et al (1991) J.Immunol. Methods 138: 47
- 20. Gaines Das, R.E. et al (1993) J.Immunol.Methods 160: 147

Reagents Plessmalaan 125 1066 CX Amsterdam The Netherlands Phone: +31 20 512 3590 Fax: +31 20 512 3570 E-mail: reagents@sanguin.nl Website: www.sanguinreagents.com





		STATIC INCUBATION	SHAKEN INCUBATION	
		Calculated mean absorbance at 450 nm		
substrate	blank	0	0	
0	pg/ml	0.014	0.019	
0.6	pg/ml	0.046	0.066	
1.9	pg/ml	0.110	0.179	
5.6	pg/ml	0.302	0.474	
16.7	pg/ml	0.793	1.245	
50	pg/ml	1.738	2.667	
150	pg/ml	2.497	> 3.000	
450	pg/ml	2.750	> 3.000	

DO NOT USE THESE DATA TO CONSTRUCT A STANDARD CURVE FOR SAMPLE VALUE CALCULATIONS