

## PeliKine compact<sup>TM</sup> human TNF $\alpha$

## Sanquin

**Specification sheet** 

REF	M1923				
Introduction	Tumour necrosis endogenous medi originally describe macrophage cyto cytotoxic factor ( TNF $\alpha$ is capable of affect the phenot that TNF $\alpha$ is not invasive stimuli ir and monocytes a response to bacte Bioassays for the been used for sex IL-1 and for this r mentioned above interference by of The Pelikine comp and specific quan cell-culture super	factor $\alpha$ (TNF $\alpha$ ) is an extreme factor of inflammatory, immune ed for their biological activities toxin (MCT), necrosin, cytoto MCF) and differentiation-indu of acting independently and ir ype and metabolism of cells i produced constitutively by no the setting of both neoplasti re thought to be the cells whi erial, viral and parasitic organi quantification of TNF $\alpha$ , inclu veral years. However, TNF $\alpha$ s eason the two commonly inte , is unaffected by IL-1, it rem ther substances. pact <sup>TM</sup> human TNF $\alpha$ ELISA kit tification of human TNF $\alpha$ in s natant.	ely potent peptide cytokine w e and host defence functions. s have been identified as TNF xin (CTX), haemorrhagic fact cing factor (DIF). a conjunction with a variety o n every tissue of the body. It rmal cells, but rather to be in c and infectious disease. In t ch contribute most to the lo sms and products. ding the cytotoxic assay on r hares many of the biological erfere in bioassays. Although ains time consuming and mig has been developed for faste erum, plasma and other body	hich serve Several s α; cachec or, macro f other fac is general duced pot his role, m cal and sy nurine fibr effects of the cytoto ht be susc er, more re fluids, as	es as an ubstances tin, phage ctors to lly thought tently by nacrophages stemic TNFα roblasts have oxic assay ceptible to eproducible s well as in
Assay procedure Kitcomponent list	See Assay proce	dure for PeliKine™ compact	ELISA kit.	T	
	Quantity	Kit component		Volume	Cap colour
	1 vial	coating antibody	100-fold concentrated	375 <i>µ</i> l	red
	1 vial	blocking reagent	50-fold concentrated	2 ml	transparent
	1 vials	TNFα standard (lyophilized)	see label	500 µl	black
	1 vial	biotinylated antibody	100-fold concentrated	375 <i>µ</i> l	yellow
	1 vial	streptavidin-poly-HRP conjugate	10,000-fold concentrated	20 <i>µ</i> l	brown
	1 bottle	HPE-dilution buffer	5-fold concentrated	55 ml	
	3 pcs	microtiter plate + lid	-	-	
	10 pcs	plate seals		-	
Sensitivity	MEAN calculated zero signal + 3 SD :1 - 3 pg/ml (shake - static incubation)2x (MEAN calculated zero signal):4 - 6 pg/ml (shake - static incubation)				
Expected values	TNF $\alpha$ values in fr	esh serum and plasma sample	es of healthy individuals are b	elow 10 p	og/ml.
Specificity	No crossreactivity was observed with the following recombinant human proteins: IL-1 $\alpha$ , IL-1 $\beta$ , IL-2, IL-3, IL-4, IL-5, IL-6, sIL-6r (GP80), IL-7, IL-8, IL-9, IL-10, IL-11, IL-13, Macrophage Colony Stimulating Factor (M-CSF), Granulocyte Colony Stimulating Factor (G-CSF), Granulocyte/Macrophage Colony Stimulating Factor (GM-CSF), Leukemia Inhibitory Factor (LIF), RANTES, Stem Cell Factor/ Mast Cell Factor (SCF/MCF), Transforming Growth Factor $\beta$ -1 (TGF $\beta$ -1), Tumour Necrosis Factor $\beta$ (TNF $\beta$ /Lymphotoxin), and Interferon $\gamma$ (IFN $\gamma$ ).				
References	<ol> <li>Aggerwal,B.B. Gutterman,J.U. (1992) Human cytokines. Blackwell Sci.Pub. ISBN 0-86542- 183-8</li> <li>Ammann,A.J. <i>et al</i> (1987) J.Clin.Immunol. <u>7</u>: 6</li> <li>Beutler,B. <i>et al</i> (1985) Science <u>229</u>: 869</li> <li>Beutler,B. <i>et al</i> (1985) Nature <u>316</u>: 552</li> <li>Carswell,E.A. <i>et al</i> (1975) Proc.Natl.Acad.Sci. <u>72</u>: 3666</li> <li>Dinarello,C.A. <i>et al</i> (1986) J.Exp.Med. <u>163</u>: 1433</li> <li>Dinarello,C.A. <i>et al</i> (1988) Rev.Inf.Diseas. <u>10</u>: 168</li> <li>Exely, A.R. <i>et al</i> (1990) Cytokine <u>2</u>: 353</li> <li>Kwiakowski,D.A. <i>et al</i> (1984) Nature <u>312</u>: 724</li> <li>Silberstein,S.B. <i>et al</i> (1986) Proc.Natl.Acad.Sci. 83: 1055</li> </ol>				
	<ol> <li>Sugarman,B</li> <li>Thomson,A.</li> <li>Tori,F.M. <i>et</i></li> </ol>	J. <i>et al</i> (1985) Science <u>230</u> : W. (1991) The cytokine hand <i>al</i> (1985) Science 229: 867	943 book. Academic Press ISBN (	)-12-6896	360-7

Reagents Plessmalian 125 1066 CX Amsterdam The Netherlands Phone: + 31 20 512 3570 Fax: + 31 20 512 3570 E-mai: regente@sanquin.nl Website: www.sanquinreagents.com

Standard	A natural human TNF $\alpha$ standard has been calibrated against the WHO International Standard (TNF $\alpha$ 87/650; National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, U.K. 1 WHO Unit = 25 pg TNF $\alpha$ ).			
	The kit contains one lyophilized vial with natural $TNF\alpha$ .			
	Reconstitute the lyophilized standard by adding 500 $\mu$ l of distilled water to the vial. Incubate for 10 minutes at room temperature and mix gently. After reconstitution the standard must be kept cold (2-8°C) and stored frozen after use (<-18°C, preferably <-70°C).			
Standard curve	Label 7 tubes, one tube for each dilutions: 1000, 333, 111, 37, 12.4, 4.1 and 1.4 pg/ml. Pipette 168 $\mu$ l of working-strength dilution buffer into the tube labelled 1000 pg/ml and 150 $\mu$ l of workingstrength dilution buffer into the other tubes. Transfer 70 $\mu$ l of the TNF $\alpha$ standard (3400 pg/ml) into the first tube labelled 1000 pg/ml, mix well and transfer 75 $\mu$ l of this dilution into the second tube labelled 333 pg/ml. Repeat the serial dilutions five more times by adding 75 $\mu$ l of the previous tube of diluted standard to the 150 $\mu$ l of dilution buffer. The standard curve will contain 1000, 333, 111, 37, 12.4, 4.1, 1.4 and 0 pg/ml (dilution buffer).			
	It is recommended to prepare two separate series for each assay.			
Samples	It is recommended to dilute the test samples at least 1:2 in working-strength dilution buffer. If high levels of TNF $\alpha$ (outside the standard curve) are expected in the test samples, additional dilutions of sample i.e. 1:10 and 1:50 should also be prepared. If rheumatoid factors are expected in samples, it is recommended to add normal mouse serum (Sanquin product M1250), final concentration in the diluted sample should be 5%.			

Typical standard curve



	STATIC INCUBATION	SHAKEN INCUBATION	
	Calculated mean absorbance at 450 nm		
substrate blank	0	0	
0 pg/ml	0.012	0.029	
1.4 pg/ml	0.016	0.046	
4.1 pg/ml	0.028	0.064	
12.4 pg/ml	0.061	0.132	
37 pg/ml	0.144	0.391	
111 pg/ml	0.403	0.895	
333 pg/ml	1.076	2.361	
1000 pg/ml	2.157	≥ 3.000	

DO NOT USE THESE DATA TO CONSTRUCT A STANDARD CURVE FOR