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YK131 Mouse/Rat CRF-HS ELISA

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FOR RESEARCH LABORATORY USE ONLY

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VIII. 文 献

1. Vale W, Spiess J, Rivier C, Rivier J: (1981) Characterization of a 41-residue ovine hypothalamic peptide that stimulates secretion of corticotropin and beta-endorphin. *Science*, **213**, 1394-1397.
2. Shibahara S, Morimoto Y, Furutani Y, Notake M, Takahashi H, Shimizu S, Horikawa S, Numa S: (1983) Isolation and sequence analysis of the human corticotropin-releasing factor precursor gene. *EMBO J*, **2**, 775-779.
3. Rivier J, Spiess J, Vale W: (1983) Characterization of rat hypothalamic corticotropin-releasing factor. *Proc Natl Acad Sci USA*, **80**, 4851-4855.
4. Seasholtz AF, Bourbonnais FJ, Harnden CE, and Camper SA: (1991) Nucleotide Sequence and Expression of the Mouse Corticotropin-Releasing Hormone Gene. *Mol. Cell. Neurosci*, **2**, 266-273
5. Vaughan J, Donaldson C, Bittencourt J, Perrin MH, Lewis K, Sutton S, Chan R, Turnbull AV, Lovejoy D, Rivier C, et al: (1995) Urocortin, a mammalian neuropeptide related to fish urotensin I and to corticotropin-releasing factor. *Nature*, **378**, 287-292.
6. Donaldson CJ, Sutton SW, Perrin MH, Corrigan AZ, Lewis KA, Rivier JE, Vaughan JM, Vale WW: (1996) Cloning and characterization of human urocortin. *Endocrinology* **137**:2167-2170
7. Chappell PB, Smith MA, Kilts CD, Bissette G, Ritchie J, Anderson C, Nemeroff CB: (1986) Alterations in corticotropin-releasing factor-like immunoreactivity in discrete rat brain regions after acute and chronic stress. *J Neurosci*, **6**, 2908-2914.
8. Nakane T, Audhya T, Hollander CS, Schlesinger DH, Kardos P, Brown C, Passarelli J: (1986) Corticotrophin-releasing factor in extra-hypothalamic brain of the mouse: demonstration by immunoassay and immunoneutralizaion of bioassayable activity. *J Endocrino*, **111**, 143-149
9. Tsuchihashi T, Yamaguchi K, Abe K, Yanaihara N, Saito S: (1992) Production of immunoreactive corticotropin-releasing hormone in various neuroendocrine tumors. *Jpn J Clin Oncol*, **22**, 232-237.
10. Nemeroff CB, Widerlov E, Bissette G, Walleus H, Karlsson I, Eklund K, Kilts CD, Loosen PT, Vale W: (1984) Elevated concentrations of CSF corticotropin-releasing factor-like immunoreactivity in depressed patients. *Science*, **226**, 1342-1344.
11. Ur E, Grossman A: (1992) Corticotropin-releasing hormone in health and disease: an update. *Acta Endocrinol (Copenh)*, **127**, 193-199. Review.
12. Geraciotti TD Jr, Orth DN, Ekhtor NN, Blumenkopf B, Loosen PT: (1992) Serial cerebrospinal fluid corticotropin-releasing hormone concentrations in healthy and depressed humans. *J Clin Endocrinol Metab*, **74**, 1325-1330.
13. May C, Rapoport SI, Tomai TP, Chrousos GP, Gold PW: (1987) Cerebrospinal fluid concentrations of corticotropin-releasing hormone (CRH) and corticotropin (ACTH) are reduced in patients with Alzheimer's disease. *Neurology*, **37**, 535-538.

14. Davis KL, Mohs RC, Marin DB, Purohit DP, Perl DP, Lantz M, Austin G, Haroutunian V: (1999) Neuropeptide abnormalities in patients with early Alzheimer disease. *Arch Gen Psychiatry*, **56**, 981-987.
15. Behan DP, Khongsaly O, Owens MJ, Chung HD, Nemeroff CB, De Souza EB: (1997) Corticotropin-releasing factor (CRF), CRF-binding protein (CRF-BP), and CRF/CRF-BP complex in Alzheimer's disease and control postmortem human brain. *J Neurochem*, **68**, 2053-2060.
16. Behan DP, Heinrichs SC, Troncoso JC, Liu XJ, Kawas CH, Ling N, De Souza EB: (1995) Displacement of corticotropin releasing factor from its binding protein as a possible treatment for Alzheimer's disease. *Nature*, **378**, 284-287.
17. Pedersen WA, McCullers D, Culmsee C, Haughey NJ, Herman JP, Mattson MP: (2001) Corticotropin-releasing hormone protects neurons against insults relevant to the pathogenesis of Alzheimer's disease. *Neurobiol D*, **8**, 492-503.
18. Delawary M, Tezuka T, Kiyama Y, Yokoyama K, Inoue T, Hattori S, Hashimoto R, Umemori H, Manabe T, Yamamoto T, Nakazawa T: (2010) NMDAR2B tyrosine phosphorylation regulates anxiety-like behavior and CRF expression in the amygdala. *Mol Brain*, **3**:37.
19. Ullrich M, Bundschu K, Benz PM, Abesser M, Freudinger R, Fischer T, Ullrich J, Renne T, Walter U, Schuh K: (2011) Identification of SPRED2 (sprouty-related protein with EVH1 domain 2) as a negative regulator of the hypothalamic-pituitary-adrenal axis. *J Biol Chem*, **286**, 9477-9488.
20. Nemoto T, Sugihara H, Mano A, Kano T, Shibasaki T: (2011) The effects of ghrelin/GHSs on AVP mRNA expression and release in cultured hypothalamic cells in rats. *Peptides*, **32**, 1281-1288.

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