



## Taq DNA Polymerase (with dNTPs), Economy

02-001 200 U  $(5\text{U/}\mu\text{l})$ , 02-001-5 5 X 200 U  $(5\text{U/}\mu\text{l})$ 

Thermus aquaticus DNA polymerase (Taq DNA polymerase) was expressed in E. coli in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and the MW is 94 kDa.. This enzyme is suitable for PCR reactions; capable of amplifying DNA with various primers.

## Applications:

1)	High-throughput	PCB
1/	TITELL UILLOUGILDUU	1 01

- 2) Colony PCR
- Incorporation of dUTP, dITP, and fluorescence-labeled nucleotides
- 4) Primer extension
- 5) Addition of a single nucleotide (adenosine) at the 3'-blunt ends

## Storage Conditions:

General composition of PCR reaction mixture (total 50µl)			
Taq DNA polymerase (5 units/µ	ıl) *0.25 μl		
10 x Reaction Buffer ( <i>Taq</i> )	5 μl		
2.5mM (each) dNTPs	$4\mu l$		
Template	<500ng		
Primer 1	$0.2 \sim 1.0 \mu M$ (final conc.)		
Primer 2	$0.2 \sim 1.0 \mu M$ (final conc.)		
Sterile distilled water	up to 50μl		
*Use of excess amount is not recommended			

20mM Tris-HCl (pH 8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.5% Tween20, 0.5% Igepal CA-630. Store at -20

Concentration:5 units/µl, where one unit is defined as the amount of enzyme that can incorporate 10 nmols of total dNTPs into an acid-insoluble material in 30 minutes at 74 when activated salmon sperm DNA was used as template/primer.

Quality Assurance: Greater than 95% purity as determined by SDS-PAGE (CBB staining) (Fig.1)

The absence of endonucleases and exonucleases was confirmed.

**PCR Test:**Good amplification result was obtained in PCR reaction using λDNA as a template (Fig.2).

## Reagents Supplied with Enzyme:

10 x Reaction Buffer (Taq): 100mM Tris-HCl (pH 8.3), 500mM KCl, 15mM MgCl<sub>2</sub>

2.5mM(each) dNTPs




Fig.1SDS-PAGE of Taq DNA polymerase

Fig.2 Amplification of DNA

Typical other

Related product: # 2-021 Pfu DNA pol#0ymerase (+dNTPs), Economy # 2-031 Pfu DNA pol#0ymerase (-dNTPs), Economy



3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013, JAPAN Tel +81(3) 5510-2347 Fax +81(3) 5510-0133

URL: http://www.sceti.jp/export/ e-mail: exp-pet@sceti.co.jp