

## ***Pfu* DNA Polymerase, Economy**

02-031 200 U (2.5U/μl), 02-031-5 5 X 200 U (2.5U/μl)

***Pyrococcus furiosus* DNA polymerase (*Pfu* DNA polymerase)** gene was expressed in *E. Coli* in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and 3' 5' exonuclease (proofreading) activity. The MW is 90 kDa, same as that of the natural *Pfu* DNA polymerase.

*Pfu* DNA polymerase is thermostable and has low error rates.

It is suitable for PCR and primer extension reactions that require high fidelity synthesis.

*Pfu* DNA polymerase-generated PCR fragments are blunt-ended.

### **Applications:**

- 1) cloning
- 2) DNA expression
- 3) site-directed mutagenesis

### **Storage Conditions:**

50mM Tris-HCl (pH 8.2), 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.1% Tween20, 0.1% Igepal CA-630

Store at -20

**Concentration:** 2.5 units/μl, where one unit is defined as the amount of enzyme that can incorporate 10 nmols of dNTPs into an acid-insoluble material in 30 minutes at 72 when activated salmon sperm DNA was used as template/primer.


**Quality Assurance:** Greater than 95% of protein determined by SDS-PAGE (CBB staining)(Fig.1)

The absence of endonucleases and exonucleases was confirmed.

**PCR Test:** Good amplification result was obtained in PCR reaction using λDNA as a template (Fig.2).

### **Reagents Supplied with Enzyme:**

10 x Reaction Buffer (*Pfu*): 200mM Tris-HCl (pH 8.8), 100mM KCl, 100mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 20mM MgSO<sub>4</sub>, 1% TritonX-100, 1 mg/ml BSA



*Pfu* DNA polymerase

### **PCR condition**

98	10sec	} 30cycles
55	30sec	
72	10min	
(2min in the case of 2kb DNA)		

lane

M : marker  
1 : 2 kb  
2 : 4 kb  
3 : 6 kb  
4 : 8 kb




Fig.1 SDS-PAGE of *Pfu* DNA polymerase

Fig.2 Amplification of λDNA

**Related products: # 02-001 Taq DNA Polymerase (+dNTPa) #02-011 Taq DNA Polymerase**