

Streptolysin O (Hemolytic streptococcus)

01-531 20 u g

Streptolysin O (SLO) is a membrane-damaging extracellular toxin produced by hemolytic streptococci. The membrane-damaging activity is measured by hemolysis of red-blood cell. SLO is oxygen-sensitive and is easily inactivated in its presence but can be reactivated by thiol compounds, so it is also called thiol-activated cytolysin (2). SLO is produced not only by Group A hemolytic streptococci but also by Group C and Group G strains. The amino acid sequences are highly conserved among them and their homology is over 98%.

The product was highly purified from E.coli over-expressing SLO of Group C hemolytic streptococci. The specific activity is as high as 1,900,000 hemolytic units (HU) /mg and forms a big hole on the cell membrane, which enables the introduction of protein inside the cells in vivo (Ref.1).

Applications

- 1) Antigen for the measurement of anti-streptolysin 0 antibody (AS0) (diagnostic reagent).
- 2) Reagent for membrane pore formation to introduce small to macromolecules into living cells (Ref.1) or deplete intracellular ATP (Ref.3).

Specifications

Measurement of the activity: Difinition of 1HU is the activation of 50% hemolysis on 3% sheep red blood cell

Purity: Over 98% is SLO by SDS-PAGE

Form: 1 mg/ml in PBS (-), 1 mM DTT, 50% glycerol, filtered-sterile Storage: -20°C (long period, -70°C). Inactivated SLO can be reactivated by thiol reagents such as 20 mM cysteine or 10 mM DTT (Ref.2)

References: This product has been used in Ref. 3

- 1. Walev I. et al. Delivery of proteins into living cells by reversible membrane permeabilization with streptolysin-O,PNAS, 98:3185 (2001)
- 2. M. Palmer: The family of thiol-activated, cholesterol-binding cytolysins. *Toxicon*, 39, 1681-1689 (2001).
- 3. Maeda Y. et al. GPHR is a novel anion channel critical for acidification and functions of the Golgi apparatus. Nature Cell Biol. 10:1135-45(2008)




Figure. SLO analyzed by SDS-PAGE.

<Distributed by >: SCETI K.K. DF Kasumigaseki Place, 3-6-7 Kasumigaseki,

Chiyoda-ku Tokyo 100-0013 JAPAN

Tel: +81-3-5510-2347 Fax: +81-3-5510-0134

E-mail: exp-pet@sceti.co.jp URL: www.sceti.co.jp/export/ <Manufactured by > : BioAcademia,Inc. 7-7-18 Saito-Asagi, Ibaraki, Osaka 567-0085, JAPAN

^{*}This product is for research use only, not for human use. It should be handled carefully to avoid injection (mouse LD50, 8 ug/kg)