



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com

Cat. No. Y-1055

EIA Buffer Concentrate

Reconstitution Information:

Each bottle of the EIA Buffer Concentrate contains 50 ml of 20x concentrated EIA buffer. *Dilute the EIA Buffer Concentrate to 1,000 ml (1 liter) with sterile deionized or USP pure water (18 MOhm) and mix well.*

Composition

The Buffer Concentrate contains Trizma base and Trizma hydrochloride (pH 7.4), NaCl, Tween 20, 0.2% Thimerosal and 1 % BSA.

Recommended Use

The EIA Buffer is used for several steps in our EIA (Enzyme Immunoassay) protocols.

Storage

Please refrigerate (0 - 4° C).