

Anti human EZFIT / ZNF71 mouse monoclonal antibody

EZFIT: Zinc finger protein 71

Code No	PP-K9716-00
0000110	
	old No. 2ZK9716H
Clone No.	K9716
Lot.	A-2
Concentration	1 mg/mL
Volume	100 uL
Ig Class	G2a
Description	

Application / Recommended Concentration

In order to obtain the best results, optimal working dilutions should be determined by each individual user.

Western Blot	2 ug/mL
Non reducing Western Blot	Not yet tested
ELISA	0.1 ug/mL
Immunoprecipitation	Not yet tested
Supershift Assay	Not yet tested
Chromatin immunoprecipitation	Not yet tested
Immunohistochemistry	Not yet tested

Nomenclature	-		
Genbank	AF269249	•	
Origin	Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human EZFIT (6-199 aa) .	Storage	Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.
Specificity	This antibody specifically recognizes human EZFIT. Not yet tested in other species.		
	,	Reference	Mataki C, et al., J Atheroscler Thromb., 7(2): 97-103, 2000.
Purification	Ammonium sulfate fractionation	-	
Formulation	Physiological saline with 0.1% NaN3 as a preservative.	Notes	Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited.

MADE IN JAPAN

Apr 23, 2008

Distributed by



3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013, Japan
Tel +81(3) 5510-2347 Fax +81(3) 5510-0133
URL: http://www.sceti.co.jp/export/ e-mail: exp-pet@sceti.co.jp



Manufactured by