

## Peninsula Laboratories, LLC

## A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071 www.bachem.com

Cat. No. S-4000 Angiotensin II

## **Specifications**

Number of Slides: 150

Chromogens: DAB (brown) and AEC (red)

Sample: Paraffin-embedded or frozen sections
Shelf Life: One year from production date

Storage: 4°C

Assay Time: 2-3 hours following primary incubation

*CROSS REACTION WITH:	%:	Kit Components
Angiotensin II	100	•30% Hydrogen peroxide Solution
[Val <sup>5</sup> ]-Angiotensin II	100	Blocking Solution
Angiotensin III	100	<ul><li>Primary Antibody</li></ul>
Renin Substrate	0.9	<ul> <li>Biotinylated Secondary Antibody</li> </ul>
Angiotensin I	0.7	<ul> <li>SA-HRP Conjugate Solution</li> </ul>
ANF (1-28) (human)	0.0	<ul> <li>DAB Concentrated Solution</li> </ul>
[Arg <sup>8</sup> ]-Vasopressin	0	•Substrate Buffer for DAB (ready to use)
[Sar <sup>1</sup> , Ile <sup>8</sup> ]-Angiotensin II	0	AEC Concentrated Solution
		•Concentrated Substrate Buffer for AEC
		<ul> <li>Hydrogen Peroxide Concentrate</li> </ul>
		Hematoxylin Counterstain
		Mounting Solution
*tastad by RIA		$\epsilon$

## \*tested by RIA