



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com

Immunohistochemistry Staining Detection Kit S-4030

Specifications

Number of Slides:	150
Chromogens:	DAB (brown), Ni-intensified DAB (dark brown), and AEC (red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation
Note::	Primary antibody sold separately

Kit Components

- 30% Hydrogen peroxide Solution
- Blocking Solution
- Biotinylated Secondary Antibody
- SA-HRP Conjugate Solution
- DAB Concentrated Solution
- Nickel-intensified DAB Concentrated Solution
- AEC Concentrated Solution
- Substrate Concentrated Buffer
- Hydrogen Peroxide Concentrate
- Hematoxylin Counterstain
- Mounting Solution