

Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071 www.bachem.com

Cat. No. \$-4032 Humanin (human)

Specifications

Number of Slides: 150

Chromogens: DAB (brown), Ni-intensified DAB (dark brown), and AEC

(red)

Sample: Paraffin-embedded or frozen sections

Shelf Life: One year from production date

Storage: 4°C

Assay Time: 2-3 hours following primary incubation

Kit Components *CROSS REACTION WITH: %: •30% Hydrogen peroxide Solution Humanin (Human) 100 •Blocking Solution β- Amyloid (1-43) (human) 0 β -Amyloid (1-42) (human) 0 Primary Antibody β-Amyloid (1-40) (human) 0 •Biotinylated Secondary Antibody •SA-HRP Conjugate Solution *tested by RIA •DAB Concentrated Solution Nickel-intensified DAB Concentrated Solution AEC Concentrated Solution •Substrate Concentrated Buffer •Hydrogen Peroxide Concentrate •Hematoxylin Counterstain Mounting Solution