

Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction	
Catalog #	RK-034-01	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	10-1280 pg/tube in 100µl Standard Solution = 100-12800 pg/ml	
Lowest Detection Limit	1103 pg/ml	
Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	Inhibin, α-Subunit (1-32) (Human)	100
	Inhibin, α-Subunit (1-32) (Porcine)	100
	Seminal Plasma Inhibin-Like Peptide (Human)	0
	Substance P	0
	ACTH (Human)	0
	Somatostatin	0
Sequence	Ser - Thr - Pro - Leu - Met - Ser - Trp - Pro - Trp - Ser - Pro - Ser - Ala - Leu - Arg - Leu - Leu - Gln - Arg - Pro - Pro - Glu - Glu - Pro - Ala - Ala - His - Ala - Asn - Cys - His - Arg	