


Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction	
Catalog #	RK-012-23	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	100-12800 pg/ml	
Lowest Detection Limit	270 pg/ml	
Cross Reactivity	Peptide	% Cross-reactivity
	Prepro-Augurin (108-132) (H)	100
	Prepro-Augurin (71-107) (H)	0
	Prepro-Augurin (133-148) (H)	0
Sequence	Asp - Arg - Asn - Gly - His - Glu - Tyr - Tyr - Gly - Asp - Tyr - Tyr - Gln - Arg - His - Tyr - Asp - Glu - Asp - Ser - Ala - Ile - Gly - Pro - Arg	

Standard Curve Image



Linear Range: 270-5324 pg/ml