


Catalog #	RK-012-25	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	10-1280 pg/ml	
Lowest Detection Limit	200 pg/ml	
Cross Reactivity	<i>Peptide</i>	<i>% Cross-reactivity</i>
	Augurin (Human) / Prepro-Augurin (71-148)	100
	Prepro-Augurin (133-148)	100
	Augurin (Human) Short Form / Prepro-Augurin (71-132)	0
	ACTH (Human)	0
	ACTH (Rat, Mouse)	0
	CRF (Human, Rat, Mouse, Canine, Feline)	0
	Vasopressin (AVP) [Arg8] (Human, Rat, Mouse, Ovine)	0
Sequence	pGlu - Leu - Trp - Asp - Arg - Thr - Arg - Pro - Glu - Val - Gln - Gln - Trp - Tyr - Gln - Gln - Phe - Leu - Tyr - Met - Gly - Phe - Asp - Glu - Ala - Lys - Phe - Glu - Asp - Asp - Ile - Thr - Tyr - Trp - Leu - Asn - Arg - Asp - Arg - Asn - Gly - His - Glu - Tyr - Tyr - Gly - Asp - Tyr - Tyr - Gln - Arg - His - Tyr - Asp - Glu - Asp - Ser - Ala - Ile - Gly - Pro - Arg - Ser - Pro - Tyr - Gly - Phe - Arg - His - Gly - Ala - Ser - Val - Asn - Tyr - Asp - Asp - Tyr	

Standard Curve
Image



Linear Range: 200-800 pg/ml