

Catalog # RK-013-01

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube in 100µl of Standard Solution =
100-12800 pg/ml

**Lowest
Detection Limit** 530.75 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
-------------------------	-----------------------	---------------------------------

	Dendroaspis Natriuretic Peptide (DNP)	100
--	--	-----

	α-ANP (1-28) (Human, Canine)	0
--	------------------------------	---

	Prepro-ANP (26-55) (Human)	0
--	----------------------------	---

	CNP-22 (Human, Porcine, Rat)	0
--	------------------------------	---

	BNP-32 (Human)	0
--	----------------	---

Sequence Glu - Val - Lys - Tyr - Asp - Pro - Cys - Phe - Gly - His - Lys
- Ile - Asp - Arg - Ile - Asn - His - Val - Ser - Asn - Leu - Gly
- Cys - Pro - Ser - Leu - Arg - Asp - Pro - Arg - Pro - Asn -
Ala - Pro - Ser - Thr - Ser - Ala ~ Disulfide bridge between
Cys7 & Cys23