

Catalog # RK-029-30

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 1-128 pg/tube of 100ul Standard Solution = 10-1280 pg/ml

Lowest Detection Limit 90.2 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	POMC Precursor (27-52) (Porcine)	100

Sequence Trp - Cys - Leu - Glu - Ser - Ser - Gln - Cys - Gln - Asp - Leu - Ser - Thr - Glu - Ser - Asn - Leu - Leu - Ala - Cys - Ile - Arg - Ala - Cys - Lys - Pro

Disulfide Bridge: Cys2 - Cys24 and Cys8 - Cys20