

Catalog #	RK-070-30
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	10-1280pg/tube in 100µl of Standard Solution = 100-12,800 pg/ml
Lowest Detection Limit	350 pg/ml
Sequence	Glu - Ile - Gly - Asp - Glu - Glu - Ser - Ala - Lys - Phe - Pro - Ile - Gly - Arg - Arg - Asp - Phe - Asp - Met - Leu - Arg - Cys - Met - Leu - Gly - Arg - Val - Tyr - Arg - Pro - Cys - Trp - Gln - Val – Disulfide bridge between Cys22 and Cys31