


Anti Human Urate Transporter 1 (URAT1) Polyclonal Antibody

Package Size	25 μ g (250 μ L / vial)
Format	Rabbit polyclonal antibody 0.1mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20°C Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with partial C-terminal peptide of human URAT1 by peptide affinity chromatography.
Working dilution	For Immunohistochemistry ; 1~5 μ g/mL For Western blotting ; 0.1~0.5 μ g/mL
HGNC* Name	SLC22A12 (solute carrier superfamily 22A12)

*HGNC: Human Gene Nomenclature Committee



Immunohistochemistry

Sample: Human proximal tubule

Preparation of antibodies and instruction:

Kanai Y, Anzai N

Department of Pharmacology and Toxicology,

Kyorin University School of Medicine

Anti Human Urate Transporter 1 (URAT1) Polyclonal Antibody

【Reference】

1. Anzai N., Miyazaki H., Noshiro R., Khamdang S., Chairoungdua A., Shin HJ., Enomoto A., Sakmoto S., Hirata T., Tomita K., Kanai Y., Endou H. : The Multivalent PDZ Domain-Containing Protein PDZK1 Regulates Transport Activity of Renal Urate-Anion Exchanger URAT1 via Its C-terminal *J Biol Chem.*2004 Aug 10 in press
2. Ichida K., Hosoyamada M., Hisatome I., Enomoto A., Hikita M., Endou H., Hosoya T. : Clinical and molecular analysis of patients with renal hypouricemia in Japan-influence of URAT1 gene on urinary urate excretion. *J Am Soc Nephrol.*Jan;15(1).164-73, 2004
3. Koepsell H., Endou H. : The SLC22 drug transporter family. *Pflugers Arch.*Feb;447(5).666-76, 2004. Epub 2003 Jul 19.
4. Anzai N., Miyazaki H., Endou H. : Molecular biology of urate transporter. *Seikagaku* Feb;76(2).101-10, 2004 Review. Japanese. No abstract available
5. Tanaka M., Itoh K., Matsushita K., Wakita N., Adachi M., Nonoguchi H., Kitamura K., Hosoyamada M., Endou H., Tomita K. : Two male siblings with hereditary renal hypouricemia and exercise-induced ARF. *Am J Kidney Dis.*Dec;42(6).1287-92, 2003
6. Enomoto A., Niwa T., Kanai Y., Endou H. : Urate transporter and renal hypouricemia. *Rinsyo Byori.* Sep;51(9).892-7.Japanese, 2003
7. Enomoto A., Kimura H., Chairoungdua A., shigeta Y., Jutabha P., Cha SH., Hosoyamada M., Takeda M., Skine T., Igarashi T., Mtsuo H., Kikuchi Y., oda T., Ichida K., Hosoya T., Shimokata K., Niwa T., Kanai Y., Endou H. : Molecular identification of a renal urate anion exchanger that regulates blood urate levels. *Nature.*May 23;417(6887).447-52, 2002. Epub 2002 Apr 14.

Supplier

SCETI
SCETI K.K.

3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013, JAPAN
Tel +81(3) 5510-2347 Fax +81(3) 5510-0133
URL: <http://www.sceti.jp/export/> e-mail: exp-pet@sceti.co.jp