

E-mail: isao.hirose@sceti.co.jp

KX447	Anti Human GPR149 Monoclonal Antibody (Clone No. 1D2)			
Primary Source	HGNC: 23627	Gene ID	344758	
Туре	Monoclonal		Keyword	
Immunogen	Partial peptide of human GPR149	G protein-coup	G protein-coupled receptor 149; IEDA; PGR10	
Raised in	GANP mouse			
Myeloma	P3U1			
Clone number	1D2			
Isotype	lgG2b,ĸ			
Source	Serum-free medium		Application	
Purification notes	ProteinG	WB	Not tested	
Cross Reactivity	Rat	IHC	5.0 μg/mL	
Concentration	0.25 mg/mL	ICC	Not tested	
Contents (Volume)	50 μg (200 μL/vial)	ELISA	1.0 μg/mL	
Label	Unlabeled	FCM	Not tested	
Buffer	PBS [containing 2 % Block Ace as a stabilizer, 0.1 % Proclin as a bacteriostat]	Neutralization	Not tested	
Storage	Store below –20 °C. Once thawed, store at 4 °C. Repeated freeze-thaw cycles should be avoided.	IP	Not tested	





This product is generated from GANP®

WARNING AND PRECAUTION

- Not for diagnostic use. The safety and efficacy of product in diagnostic or other clinical uses has not been catabilished.
- 2. Harmful by inhalation, in contact with skin and if swallowed. Do not breathe dust. Avoid contact with skin and eyes.
- If contact with skin and eyes, wash all affected areas with large volume of water. If inhaled remove to fresh air. In severe case obtain medical attention.
- 4. Wash hand thoroughly after handling the product.

SCETI

- 5. Do not use this product if container is broken or some contaminants are detected.
- 6. When preserving the product, Close the container, ensure it does not fall aside or down.
- 7. Dispose of the container and expired reagents in accordance with federal, state and local government regulations.
- 8. Do not use the container and accessories of the product for other purpose.

取り扱い上の注意

この添付文書をよく読んでから使用して下さい。

- 1. 本品は研究用試薬であり、医薬品その他の目的にはご使用になれません。
- 2. 取り扱い中は皮膚、粘膜、着衣に触れたり、目に入らないように適切な措置を行って下さい。
- 3. 試薬が誤って目やロに入った場合には、水で十分に洗い流すなどの応急処置を行い、必要があれば医師の手当を受けて下さい。
- 4. 取り扱い後には手洗いを十分に行って下さい。
- 5. 容器の破損、異物混入等異常が認められた物は使用しないで下さい。
- 6. 試薬を保管する場合は、蓋をし、転倒落下防止を確実にし、指定の貯蔵方法で保管して下さい。
- 7. 使用後の容器は、廃棄物に関する規定に従って処理して下さい。
- 8. 容器、付属品等の他目的への転用は保証できません。