



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com



Cat. No. S-1153 EIA for Substance P (for human serum and plasma samples)

SPECIFICATIONS: Extraction-free kit for human serum/plasma


STANDARD: 1 μ g H-1890 H-Arg-Pro-Lys-Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-NH₂

PROTOCOL: V (Ab1hr.Std2hr.BtON)

RANGE: 0-25 ng/ml

AVERAGE IC₅₀: 0.5 ng/ml

SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.02 to 25 ng/ml
Stock	1000	
S1	25.00	Add 25 μ l stock + 975 μ l diluent
S2	6.25	Add 200 μ l S1 + 600 μ l diluent
S3	1.56	Add 200 μ l S2 + 600 μ l diluent
S4	0.39	Add 200 μ l S3 + 600 μ l diluent
S5	0.10	Add 200 μ l S4 + 600 μ l diluent
S6	0.02	Add 200 μ l S5 + 600 μ l diluent
S0	0.00	500 μ l diluent

CROSS REACTION:

PEPTIDE	%
Substance P	100
Substance P (2-11)	100
Substance P (3-11)	100
Substance P (4-11)	100
Substance P (5-11)	100
Substance P (6-11)	<5
Substance P (7-11)	<0.01
Neuropeptide K (porcine)	<0.01
Neurokinin A	<0.01
Neurokinin B	0
Endothelin-1 (human, bovine, canine, mouse, porcine, rat)	<0.5
Physalaemin	0