



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com



Cat. No. S-1157

EIA for Endothelin-1 (human, bovine, dog, mouse, porcine, rat) (For rat Serum and Plasma Samples)

SPECIFICATIONS: Extraction-free kit for rat serum/plasma


STANDARD: 1 ug H-6995 H-Cys-Ser-Cys-Ser-Ser-Leu-Met-Asp-Lys-Glu-Cys-Val-Tyr-Phe-Cys-His-Leu-Asp-Ile-Ile-Trp-OH, (Disulfide bonds between Cys1 and Cys15/Cys3 and Cys11)

PROTOCOL: V (Ab1hr.Std2hr.BtON)

RANGE: 0-25 ng/ml

AVERAGE IC₅₀: 0.4 ng/ml

SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.02 to 25 ng/ml
Stock	1000	
S1	25.00	Add 25 µl stock + 975 µl diluent
S2	6.25	Add 200 µl S1 + 600 µl diluent
S3	1.56	Add 200 µl S2 + 600 µl diluent
S4	0.39	Add 200 µl S3 + 600 µl diluent
S5	0.10	Add 200 µl S4 + 600 µl diluent
S6	0.02	Add 200 µl S5 + 600 µl diluent
S0	0.00	500 µl diluent

CROSS REACTION:

PEPTIDE	%
Endothelin-1 (human, porcine, canine, rat, mouse, bovine)	100
Endothelin-2 (human, canine)	100
Big Endothelin-1 (human)	100
Endothelin-3 (human, rat, porcine, rabbit)	0.005
Big Endothelin-1 (porcine)	100
AANF (1-28) (human)	<0.01
BNP-32 (human)	0
Angiotensin II	0
(Arg ⁸)-Vasopressin	0
VIP (human, bovine, porcine, rat)	0