



Peninsula Laboratories, LLC
A Member of the Bachem Group
 305 Old County Road, San Carlos, CA 94070
 Tel: (800) 922-1516 • (650) 592-5392
 Fax: (650) 595-4071
 www.bachem.com




Cat. No. S-1204
Dynorphin-A (1-8)

SPECIFICATIONS: Formulated for extracted sample

STANDARD: 1 µg H-2655 H-Tyr-Gly-Gly-Phe-Leu-Arg-Arg-Ile-OH

PROTOCOL: III (Std.Ab1hr.Bt)
RANGE: 0-5 ng/ml
AVERAGE IC₅₀: 0.13 ng/ml
SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.00 to 5 ng/ml
Stock	1000	
S1	5.00	Add 5 µl stock + 995 µl diluent
S2	1.25	Add 200 µl S1 + 600 µl diluent
S3	0.31	Add 200 µl S2 + 600 µl diluent
S4	0.08	Add 200 µl S3 + 600 µl diluent
S5	0.02	Add 200 µl S4 + 600 µl diluent
S6	0.00	Add 200 µl S5 + 600 µl diluent
S0	0.00	500 µl diluent

CROSS REACTION:

PEPTIDE	%
Dynorphin A (1-8)	100
Dynorphin A	0.01
Prodynorphin (209-240) (porcine)	0.01
Dynorphin A (1-13)	0.01
Dynorphin B	0
α-Neoendorphin	0
β-Endorphin (human)	0
Leu-Enkephalin	0
Met-Enkephalin	0