



Peninsula Laboratories, LLC
A Member of the Bachem Group
 305 Old County Road, San Carlos, CA 94070
 Tel: (800) 922-1516 • (650) 592-5392
 Fax: (650) 595-4071
 www.bachem.com




Cat. No. S-1215
Inhibin α -Subunit (1-32) (porcine)

SPECIFICATIONS: Formulated for extracted samples

STANDARD: 1 μ g (not available separately) H-Ser-Thr-Ala-Pro-Leu-Pro-Trp-Pro-Trp-Ser-Pro-Ala-Ala-Leu-Arg-Leu-Leu-Gln-Arg-Pro-Pro-Glu-Glu-Pro-Ala-Val-His-Ala-Asp-Cys-His-Arg-OH

PROTOCOL: III (Std.Ab1hr.Bt)
RANGE: 0-5 ng/ml
AVERAGE IC₅₀: 0.12 ng/ml
SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.00 to 5 ng/ml
Stock	1000	
S1	5.00	Add 5 μ l stock + 995 μ l diluent
S2	1.25	Add 200 μ l S1 + 600 μ l diluent
S3	0.31	Add 200 μ l S2 + 600 μ l diluent
S4	0.08	Add 200 μ l S3 + 600 μ l diluent
S5	0.02	Add 200 μ l S4 + 600 μ l diluent
S6	0.00	Add 200 μ l S5 + 600 μ l diluent
S0	0.00	500 μ l diluent

CROSS REACTION:

PEPTIDE	%
Inhibin, α -Subunit 1-32 (porcine)	100
Inhibin, α -Subunit 1-32 (human)	100
α -Inhibin-31 (human)	0
Substance P	0
ACTH (1-39) (human)	0
Somatostatin-14	0