



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com



Cat. No. S-1279

EIA for Neuromedin S (human) (for human serum and plasma samples)

SPECIFICATIONS: Extraction-free kit for human serum/plasma


STANDARD: 1 µg H-6232 H-Ile-Leu-Gln-Arg-Gly-Ser-Gly-Thr-Ala-Ala-Val-Asp-Phe-Thr-Lys-Lys-Asp-His-Thr-Ala-Thr-Trp-Gly-Arg-Pro-Phe-Phe-Leu-Phe-Arg-Pro-Arg-Asn-NH₂

PROTOCOL: V (Ab1hr.Std2hr.BtON)

RANGE: 0-10 ng/ml

AVERAGE IC₅₀: 0.6 ng/ml

SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.01 to 10 ng/ml
Stock	1000	
S1	10.00	Add 10 µl stock + 990 µl diluent
S2	2.50	Add 200 µl S1 + 600 µl diluent
S3	0.63	Add 200 µl S2 + 600 µl diluent
S4	0.16	Add 200 µl S3 + 600 µl diluent
S5	0.04	Add 200 µl S4 + 600 µl diluent
S6	0.01	Add 200 µl S5 + 600 µl diluent
S0	0.00	500 µl diluent

CROSS REACTION:

PEPTIDE	%
Neuromedin S (human)	100
Neuromedin S (rat)	25
Neuromedin U-25 (human)	0
Neuromedin U-23 (rat)	0
PrePro-Neuromedin U (104-136) (human)	0
PrePro-Neuromedin S (70-103) (human)	0
Neuropeptide Y (human, rat)	0
Ghrelin (human)	0