



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com




Cat. No. S-1283 D-Marker peptide

SPECIFICATIONS: Formulated for extracted samples

STANDARD: 1 µg (not available separately) H-D-Cys-Gly-D-Pro-Gly-Gly-D-Tyr-D-Lys-D-Val-D-Leu-D-Val-D-Leu-D-Asn-D-Pro-D-Ser-D-Val-D-Ala-D-Ala-D-Thr-OH

PROTOCOL: 1 (Std.Ab.Bt)
RANGE: 0-10 ng/ml
AVERAGE IC₅₀: 0.4 ng/ml
SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.01 to 10 ng/ml
Stock	1000	
S1	10.00	Add 10 µl stock + 990 µl diluent
S2	2.50	Add 200 µl S1 + 600 µl diluent
S3	0.63	Add 200 µl S2 + 600 µl diluent
S4	0.16	Add 200 µl S3 + 600 µl diluent
S5	0.04	Add 200 µl S4 + 600 µl diluent
S6	0.01	Add 200 µl S5 + 600 µl diluent
S0	0.00	500 µl diluent

CROSS REACTION:

PEPTIDE

D-Marker peptide
L-Marker peptide (CGPGGYKVLVLNPSVAAT)

%

100
0