



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com



Cat. No. S-1338

Calcitonin (rat)

(For rat serum and plasma samples)

SPECIFICATIONS: Extraction-free kit for rat serum/plasma


STANDARD: 1 µg H-3072 H-Cys-Gly-Asn-Leu-Ser-Thr-Cys-Met-Leu-Gly-Thr-Tyr-Thr-Gln-Asp-Leu-Asn-Lys-Phe-His-Thr-Phe-Pro-Gln-Thr-Ser-Ile-Gly-Val-Gly-Ala-Pro-NH₂, (Disulfide bond)

PROTOCOL: V (Ab1hr.Std2hr.BtON)

RANGE: 0-25 ng/ml

AVERAGE IC₅₀: 0.65 ng/ml

SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.02 to 25 ng/ml
Stock	1000	
S1	25.00	Add 25 µl stock + 975 µl diluent
S2	6.25	Add 200 µl S1 + 600 µl diluent
S3	1.56	Add 200 µl S2 + 600 µl diluent
S4	0.39	Add 200 µl S3 + 600 µl diluent
S5	0.10	Add 200 µl S4 + 600 µl diluent
S6	0.02	Add 200 µl S5 + 600 µl diluent
S0	0.00	500 µl diluent

CROSS REACTION:

PEPTIDE	%
Calcitonin (rat)	100
Calcitonin (human)	87
Calcitonin (salmon I)	0
Calcitonin (eel)	0
Calcitonin (chicken)	0
α-CGRP (human)	0.1
α-CGRP (rat)	0
Amylin (rat)	0
Amylin (free acid) (human)	0
Calcitonin C-Terminal Adjacent Peptide (rat)	0
Calcitonin C-Terminal Flanking Peptide (human)	0

Cross-Reactivity data is taken from Enzyme immunoassay applications.