



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com



Cat. No. S-1363

(Lys⁸)-Vasopressin (Terlipressin (4-12)) (For rat serum and plasma samples)

SPECIFICATIONS: Extraction-free kit for rat serum/plasma


STANDARD: 1 µg H-2530 H-Cys-Tyr-Phe-Gln-Asn-Cys-Pro-Lys-Gly-NH₂ (Disulfide bond)

PROTOCOL: V (Ab 1hr.Std 2hr.BtON)

RANGE: 0-100 ng/ml

AVERAGE IC₅₀: 1.5 ng/ml

SHELF-LIFE: One year from production date



Suggested Preparation of Standards		
	ng/ml	Range: 0.10 to 100 ng/ml
Stock	1000	
S1	100.00	Add 100 µl stock + 900 µl diluent
S2	25.00	Add 200 µl S1 + 600 µl diluent
S3	6.25	Add 200 µl S2 + 600 µl diluent
S4	1.56	Add 200 µl S3 + 600 µl diluent
S5	0.39	Add 200 µl S4 + 600 µl diluent
S6	0.10	Add 200 µl S5 + 600 µl diluent
S0	0.00	500 µl diluent

CROSS REACTION:

PEPTIDE

%

(Lys⁸)-Vasopressin

100

(Arg⁸)-Vasopressin

100