



## PENINSULA LABORATORIES, LLC

**A Member of the Bachem Group**

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

[www.bachem.com](http://www.bachem.com)

**Cat. No. S-3062**

### **Galanin (porcine)**

#### **Specifications**

Number of Slides: 150  
Sample: Paraffin-embedded or frozen sections  
Shelf Life: One year from production date  
Storage: 4°C  
Assay Time: 2-3 hours following primary incubation

#### **Kit Components**

- Primary antibody
- FITC-labeled secondary antibody
- Normal rabbit serum
- 5X phosphate buffered saline, pH 7.4 with 0.1% sodium azide
- 3% triton X-100/water
- Glycerin-PBS (9:1)
- Normal goat serum

#### **Cross Reaction:**

| <b>Peptide</b>            | <b>%</b> |
|---------------------------|----------|
| Galanin (porcine)         | 100      |
| Galanin (rat)             | 86       |
| Galanin (human)           | 0        |
| Secretin (porcine)        | 0        |
| VIP (human, porcine, rat) | 0        |
| PHM-27 (human)            | 0        |
| NPY (porcine)             | 0        |
| GMAP 1-41 amide (porcine) | 0        |
| GMAP 44-59 (porcine)      | 0        |
| Insulin (human)           | 0        |
| Substance P               | 0        |

Cross-Reactivity data is taken from Radioimmunoassay applications