



PENINSULA LABORATORIES, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com

Cat. No. S-3156

Dynorphin B (porcine)

Specifications

Number of Slides: 150
Sample: Paraffin-embedded or frozen sections
Shelf Life: One year from production date
Storage: 4°C
Assay Time: 2-3 hours following primary incubation

Kit Components

- Primary antibody
- FITC-labeled secondary antibody
- Normal rabbit serum
- 5X phosphate buffered saline, pH 7.4 with 0.1% sodium azide
- 3% triton X-100/water
- Glycerin-PBS (9:1)
- Normal goat serum

Cross Reaction:

Peptide	%
Dynorphin B (porcine)	100
Dynorphin B (1-29)	12
Dynorphin B (15-29)	0
Dynorphin A (porcine)	0
Dynorphin A (1-8) (porcine)	0
α -Neo-Endorphin (porcine)	0
β -Neo-Endorphin	0
Leu-Enkephalin	0
Met-Enkephalin	0

Cross-Reactivity data is taken from Radioimmunoassay applications