


For research use only

Anti HMGb3 Polyclonal Antibody

HMGb3 : High mobility group box 3

Package Size	25μg (100μL/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Below -20°C
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of mouse HMGb3 by peptide affinity chromatography.
Working dilution	Western Blotting : 1~3μg/mL



Western Blotting

Sample : Mouse cerebrum homogenate

Preparation of antibodies and instruction

Dr. Shigeyoshi Y. and Dr. Takashima N. at Department of Anatomy and Neurobiology, Kinki University School of Medicine

Anti HMGb3 Polyclonal Antibody

HMGb3 : High mobility group box 3

【Reference】

1. Muller S. et al. : EMBO J. 2001 Aug 15;20(16):4337-40
2. Thomas JO. et al. : Biochem Soc Trans. 2001 Aug;29(Pt 4):395-401
3. Vaccari T. et al. : Genomics. 1998 Apr 15;49(2):247-52
4. Nemeth MJ. et al. : Blood. 2003 Aug 15;102(4):1298-306
5. Ronfani L. et al. : Development. 2001 Apr;128(8):1265-73
6. Ueda HR. et al. : Nature. 2002 Aug 1;418(6897):534-9

Supplier

SCETI
SCETI K.K.

3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013, JAPAN
Tel +81(3) 5510-2347 Fax +81(3) 5510-0133
URL: <http://www.sceti.jp/export/> e-mail: exp-pet@sceti.co.jp