


Anti Mouse DC-STAMP Polyclonal Antibody


DC-STAMP(FIND) : Dendritic Cell-Specific Transmembrane Protein (IL-Four Induced)

Package Size	25μg (100μL / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20℃ Once thawed, store at 4℃. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of mouse DC-STAMP(FIND) by peptide affinity chromatography.
Working dilution	For Immunocytochemistry ; 4μg/mL

① Anti DC-STAMP antibody



② Control rabbit IgG



Immunocytochemistry

Sample : RAW-D cells stimulated with RANKL+TNF- α

①Anti DC-STAMP antibody

②Control rabbit IgG

Preparation of antibodies and instruction :

Dr. Nomiya H. at Department of Molecular Enzymology ,
Kumamoto University Graduate school of Medical Sciences

Anti Mouse DC-STAMP Polyclonal Antibody**DC-STAMP(FIND) : Dendritic Cell-Specific Transmembrane Protein (IL-Four Induced)****【Reference】**

1. Hartgers F. C. et al. : Eur. J. Immunol. 2000 Dec;30(12):3585-90.
2. Hartgers F. C. et al. : Immunogenetics. 2001 Mar;53(2):145-9.
3. Staeger H. et al. : Immunogenetics. 2001 Mar;53(2):105-13.
4. Rho J. et al.: DNA Cell Biol. 2002 Aug;21(8):541-9.
5. Kukita T, et al. : J Exp Med. 2004 Oct 4;200(7):941-6. *
6. Yagi M. et al.: J Exp Med. 2005 Aug 1;202(3):345-51.

*It is the article that used this antibody.

Supplier**SCETI**
SCETI K.K.

3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013, JAPAN

Tel +81(3) 5510-2347 Fax +81(3) 5510-0133

URL: <http://www.sceti.jp/export/> e-mail: exp-pet@sceti.co.jp