



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com



Cat. No. S-2207

Big Endothelin-1 (porcine)

SPECIFICATIONS: Formulated for extracted samples

STANDARD: 12.8 µg H-9135 H-Cys-Ser-Cys-Ser-Ser-Leu-Met-Asp-Lys-Glu-Cys-Val-Tyr-Phe-Cys-His-Leu-Asp-Ile-Ile-Trp-Val-Asn-Thr-Pro-Glu-His-Ile-Val-Pro-Tyr-Gly-Leu-Gly-Ser-Pro-Ser-Arg-Ser-OH, (Disulfide bonds between Cys1 and Cys15/Cys3 and Cys11)

TRACER: H-4964

RANGE: 0.01-1.28 ng/ml (1-128 pg/tube)

AVERAGE IC₅₀: 0.06 ng/ml (6 pg/tube)

TEST CAPACITY: 125 test tubes

ASSAY TIME: 3 days

SHELF LIFE: Six weeks after date of preparation when stored at -20°C

CROSS REACTION:

PEPTIDE:	%:
Big Endothelin-1 (porcine)	100
Endothelin-2 (human, canine)	100
Big Endothelin-1 (human)	100
Big Endothelin-1 (rat)	100
Endothelin-1 (human, porcine, canine, rat, bovine, mouse)	100
VIC (Endothelin-Beta) (mouse)	100
Endothelin-3 (human, rat, porcine, mouse rabbit)	70
Sarafotoxin B	<0.1
Atrial Natriuretic Factor (1-28) (human)	<0.01
Big Endothelin-1 fragment (22-38) (human)	0
BNP-32 (human)	0
Angiotensin I	0
Angiotensin II	0
[Arg ⁸]-Vasopressin	0
VIP (human, porcine, rat, bovine)	0
ACTH (1-39) (human)	0