



Peninsula Laboratories, LLC
A Member of the Bachem Group
305 Old County Road, San Carlos, CA 94070
Tel: (800) 922-1516 • (650) 592-5392
Fax: (650) 595-4071
www.bachem.com

Cat. No. S-4003
 β -Endorphin (rat)

Specifications

Number of Slides:	150
Chromogens:	DAB (brown), Ni-intensified DAB (dark brown), and AEC (red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation

* CROSS REACTION WITH:	%:	Kit Components
β -Endorphin (rat)	100	•30% Hydrogen peroxide Solution
β -Endorphin (porcine)	100	•Blocking Solution
β -Endorphin (equine)	100	•Primary Antibody
β -Endorphin (bovine, camel, mouse, ovine)	100	•Biotinylated Secondary Antibody
α -Endorphin	100	•SA-HRP Conjugate Solution
γ -Endorphin	60	•DAB Concentrated Solution
β -Endorphin (human)	40	•Nickel-intensified DAB Concentrated Solution
ACTH (1-39) (human)	0	•AEC Concentrated Solution
Met-Enkephalin	0	•Substrate Concentrated Buffer
α -MSH	0	•Hydrogen Peroxide Concentrate
PACAP-38 (human, mouse, ovine, porcine, rat)	0	•Hematoxylin Counterstain
		•Mounting Solution

***tested by RIA**