



Peninsula Laboratories, LLC
A Member of the Bachem Group
305 Old County Road, San Carlos, CA 94070
Tel: (800) 922-1516 • (650) 592-5392
Fax: (650) 595-4071
www.bachem.com

Cat. No. S-4008
Somatostatin-14

Specifications

Number of Slides:	150
Chromogens:	DAB (brown), Ni-intensified DAB (dark brown), and AEC (red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation

*CROSS REACTION WITH:	%:	Kit Components
Somatostatin-14	100	•30% Hydrogen peroxide Solution
Somatostatin-28	100	•Blocking Solution
[Des-Ala ¹]-Somatostatin	100	•Primary Antibody
[D-Trp ⁸]-Somatostatin	0	•Biotinylated Secondary Antibody
Prosomatostatin-32	0	•SA-HRP Conjugate Solution
Somatostatin Analog, RC-160	0	•DAB Concentrated Solution
Somatostatin Analog (CTOP-NH ₂)	0	•Nickel-intensified DAB Concentrated Solution
VIP (human, bovine, rat, porcine)	0	•AEC Concentrated Solution
NPY (porcine)	0	•Substrate Concentrated Buffer
Insulin (human)	0	•Hydrogen Peroxide Concentrate
Glucagon (1-29) (human)	0	•Hematoxylin Counterstain
Amylin (human)	0	•Mounting Solution

***tested by RIA**