



Peninsula Laboratories, LLC
A Member of the Bachem Group
305 Old County Road, San Carlos, CA 94070
Tel: (800) 922-1516 • (650) 592-5392
Fax: (650) 595-4071
www.bachem.com

Cat. No. S-4009
Substance P

Specifications

Number of Slides: 150
Chromogens: DAB (brown) and AEC (red)
Sample: Paraffin-embedded or frozen sections
Shelf Life: One year from production date
Storage: 4°C
Assay Time: 2-3 hours following primary incubation

*CROSS REACTION WITH:	%:	Kit Components
Substance P	100	•30% Hydrogen peroxide Solution
Substance P 2-11	100	•Blocking Solution
Substance P 3-11	74.5	•Primary Antibody
Substance P 4-11	12.4	•Biotinylated Secondary Antibody
Substance P 5-11	12	•SA-HRP Conjugate Solution
Substance P 6-11	4.8	•DAB Concentrated Solution
Substance P 7-11	0	•Substrate Buffer for DAB (ready to use)
Substance P 8-11	0	•AEC Concentrated Solution
Substance P 9-11	0	•Concentrated Substrate Buffer for AEC
Neurokinin A	0	•Hydrogen Peroxide Concentrate
Neurokinin B	0	•Hematoxylin Counterstain
Endothelin-1 (human, bovine, canine, mouse, porcine, rat)	0	•Mounting Solution

***tested by RIA**