



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070

Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071

www.bachem.com

Cat. No. S-4017

VIP

(human, bovine, porcine, rat)

Specifications

Number of Slides:	150
Chromogens:	DAB (brown), Ni-intensified DAB (dark brown), and AEC (red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation

*CROSS REACTION WITH:	%:	Kit Components
VIP (human, bovine, porcine, rat)	100	•30% Hydrogen peroxide Solution
VIP (guinea pig)	100	•Blocking Solution
VIP (10-28) (human, porcine, rat)	4.4	•Primary Antibody
VIP (chicken)	< 1	•Biotinylated Secondary Antibody
VIP (1-12) (human, bovine, porcine, rat)	0	•SA-HRP Conjugate Solution
PACAP-27 (human, mouse, Ovine, porcine, rat)	0	•DAB Concentrated Solution
Substance P	0	•Nickel-intensified DAB Concentrated Solution
Endothelin-1 (human, porcine, rat, canine, mouse, bovine)	0	•AEC Concentrated Solution
PHM-27 (human)	0	•Substrate Concentrated Buffer
Secretin (porcine)	0	•Hydrogen Peroxide Concentrate
Glucagon (1-29) (human)	0	•Hematoxylin Counterstain
Galanin (rat)	0	•Mounting Solution
Somatostatin-14	0	

*tested by RIA