

Peninsula Laboratories, LLC

A Member of the Bachem Group 305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392 Fax: (650) 595-4071 www.bachem.com

> Cat. No. S-4019 **Dynorphin A**

Specifications

Number of Slides:	150
Chromogens:	DAB (brown), Ni-intensified DAB (dark brown), and AEC
	(red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation

*CROSS REACTION WITH

*CROSS REACTION WITH: Dynorphin A Dynorphin A (1-13) Dynorphin A (1-8) Dynorphin B Leu-Enkephalin β-Endorphin (human) α-Neoendorphin	%: 100 42 0 0 0 0 0	Kit Components •30% Hydrogen peroxide Solution •Blocking Solution •Primary Antibody •Biotinylated Secondary Antibody •SA-HRP Conjugate Solution •DAB Concentrated Solution •Nickel-intensified DAB Concentrated Solution •AEC Concentrated Solution •Substrate Concentrated Buffer •Hydrogen Peroxide Concentrate
*tested by RIA		Hematoxylin CounterstainMounting Solution